



**QUANTUM®**

## ILLUSTRATED PARTS BREAKDOWN



**1. Cushion Assy's**

Stealth Essence Plus Complete Cushion Assy w/ 4-way Stretch Cover . . . . . 1  
Stealth Essence Foam Complete Cushion Assy w/ 4-way Stretch Cover . . . . . 3  
Stealth Simplicity G Complete Cushion Assy . . . . . 5  
Stealth Solution Complete Cushion Assy w/ 4-way Stretch Cover . . . . . 7  
Stealth Spectrum Foam Complete Cushion Assy w/ 4-way Stretch Cover . . . . . 11  
Stealth Spectrum Gel Complete Cushion Assy w/ 4-way Stretch Cover . . . . . 15  
Stealth TRU-Comfort 2 Complete Cushion Assy w/ 4-way Stretch Cover . . . . . 19

**2. Outer Covers**

Stealth Essence Foam 4-way Stretch Outer Cushion Cover . . . . . 21  
Stealth Essence Plus 4-way Stretch Outer Cushion Cover . . . . . 23  
Stealth Simplicity G Tek2 Outer Cushion Cover . . . . . 25  
Stealth Solution 4-way Stretch Outer Cushion Cover . . . . . 27  
Stealth Spectrum Foam 4-way Stretch Outer Cushion Cover . . . . . 31  
Stealth Spectrum Gel 4-way Stretch Outer Cushion Cover . . . . . 35  
Stealth TRU-Comfort 2 4-way Stretch Outer Cushion Cover . . . . . 39

**3. Inner Covers**

Urethane Inner Liner - 10"W - 19"W - Stealth Solution, Spectrum Gel, Spectrum Foam . . . . . 41  
Urethane Inner Liner - 20"W - 32"W - Stealth Solution, Spectrum Gel, Spectrum Foam . . . . . 43  
Polyurethane Inner Liner - Essence . . . . . 47  
Polyurethane Inner Liner 10"W-19"W - Stealth Solution, Spectrum Gel, & Spectrum Foam . . . . . 49  
Polyurethane Inner Liner 20"W-32"W - Stealth Solution, Spectrum Gel, & Spectrum Foam . . . . . 51  
Polyurethane Inner Liner - Stealth TRU-Comfort 2 . . . . . 55

**4. ABS Pans - Solution, Spectrum Gel, Spectrum Foam**

ABS Pan 10"W- 19"W, Stealth Simplicity, Solution, Spectrum Gel, Spectrum Foam . . . . . 57  
ABS Pan 20"W - 32"W, Stealth Simplicity, Solution, Spectrum Gel, Spectrum Foam . . . . . 59

**5. ACM Board**

ACM Boards . . . . . 61

**6. Stealth Essence Wedge Assy's**

Stealth Essence Wedge Assy's . . . . . 63



## 1. Cushion Assy's

Stealth Essence Plus Complete Cushion Assy	
Size	Ref# w/ Large Wedges
16"W x 16"D	A1
16"W x 18"D	B1
18"W x 16"D	C1
18"W x 18"D	D1
18"W x 20"D	E1
20"W x 18"D	F1
20"W x 20"D	G1

\*ABS pan non-serviceable. Included in top level assembly.



Item	Part No.	Name	Qty
A1	ST-ESSPLSPP1616L	STEALTH,CUSHION,ASSY,ESSENCE PLUS,SKIN PROTECTION AND POSITIONING,16"W,16"D,W/ L WEDGE	1
B1	ST-ESSPLSPP1618L	STEALTH,CUSHION,ASSY,ESSENCE PLUS,SKIN PROTECTION AND POSITIONING,16"W,18"D,W/ L WEDGE	1
C1	ST-ESSPLSPP1816L	STEALTH,CUSHION,ASSY,ESSENCE PLUS,SKIN PROTECTION AND POSITIONING,18"W,16"D,W/ L WEDGE	1
D1	ST-ESSPLSPP1818L	STEALTH,CUSHION,ASSY,ESSENCE PLUS,SKIN PROTECTION AND POSITIONING,18"W,18"D,W/ L WEDGE	1
E1	ST-ESSPLSPP1820L	STEALTH,CUSHION,ASSY,ESSENCE PLUS,SKIN PROTECTION AND POSITIONING,18"W,20"D,W/ L WEDGE	1
F1	ST-ESSPLSPP2018L	STEALTH,CUSHION,ASSY,ESSENCE PLUS,SKIN PROTECTION AND POSITIONING,20"W,18"D,W/ L WEDGE	1
G1	ST-ESSPLSPP2020L	STEALTH,CUSHION,ASSY,ESSENCE PLUS,SKIN PROTECTION AND POSITIONING,20"W,20"D,W/ L WEDGE	1

## 1. Cushion Assy's

Stealth Essence Foam Complete Cushion Assy			
Size	Ref# w/o Wedges	Ref# w/ Small Wedges	Ref# w/ Med Wedges
16"W x 16"D	A1	A2	A3
16"W x 18"D	B1	B2	B3
18"W x 16"D	C1	C2	C3
18"W x 18"D	D1	D2	D3
18"W x 20"D	E1	E2	E3
20"W x 18"D	F1	F2	F3
20"W x 20"D	G1	G2	G3


\*ABS pan non-serviceable. Included in top level assembly.



Item	Part No.	Name	Qty
A1	ST-ESSSPP1616	STEALTH,CUSHION,ASSY,ESSENCE,SKIN PROTECTION AND POSITIONING,16"W,16"D	1
A2	ST-ESSSPP1616S	STEALTH,CUSHION,ASSY,ESSENCE,SKIN PROTECTION AND POSITIONING,16"W,16"D,W/ S WEDGE	1
A3	ST-ESSSPP1616M	STEALTH,CUSHION,ASSY,ESSENCE,SKIN PROTECTION AND POSITIONING,16"W,16"D,W/ M WEDGE	1
B1	ST-ESSSPP1618	STEALTH,CUSHION,ASSY,ESSENCE,SKIN PROTECTION AND POSITIONING,16"W,18"D	1
B2	ST-ESSSPP1618S	STEALTH,CUSHION,ASSY,ESSENCE,SKIN PROTECTION AND POSITIONING,16"W,18"D,W/ S WEDGE	1
B3	ST-ESSSPP1618M	STEALTH,CUSHION,ASSY,ESSENCE,SKIN PROTECTION AND POSITIONING,16"W,18"D,W/ M WEDGE	1
C1	ST-ESSSPP1816	STEALTH,CUSHION,ASSY,ESSENCE,SKIN PROTECTION AND POSITIONING,18"W,16"D	1
C2	ST-ESSSPP1816S	STEALTH,CUSHION,ASSY,ESSENCE,SKIN PROTECTION AND POSITIONING,18"W,16"D,W/ S WEDGE	1
C3	ST-ESSSPP1816M	STEALTH,CUSHION,ASSY,ESSENCE,SKIN PROTECTION AND POSITIONING,18"W,16"D,W/ M WEDGE	1
D1	ST-ESSSPP1818	STEALTH,CUSHION,ASSY,ESSENCE,SKIN PROTECTION AND POSITIONING,18"W,18"D	1
D2	ST-ESSSPP1818S	STEALTH,CUSHION,ASSY,ESSENCE,SKIN PROTECTION AND POSITIONING,18"W,18"D,W/ S WEDGE	1
D3	ST-ESSSPP1818M	STEALTH,CUSHION,ASSY,ESSENCE,SKIN PROTECTION AND POSITIONING,18"W,18"D,W/ M WEDGE	1
E1	ST-ESSSPP1820	STEALTH,CUSHION,ASSY,ESSENCE,SKIN PROTECTION AND POSITIONING,18"W,20"D	1
E2	ST-ESSSPP1820S	STEALTH,CUSHION,ASSY,ESSENCE,SKIN PROTECTION AND POSITIONING,18"W,20"D,W/ S WEDGE	1
E3	ST-ESSSPP1820M	STEALTH,CUSHION,ASSY,ESSENCE,SKIN PROTECTION AND POSITIONING,18"W,20"D,W/ M WEDGE	1
F1	ST-ESSSPP2018	STEALTH,CUSHION,ASSY,ESSENCE,SKIN PROTECTION AND POSITIONING,20"W,18"D	1
F2	ST-ESSSPP2018S	STEALTH,CUSHION,ASSY,ESSENCE,SKIN PROTECTION AND POSITIONING,20"W,18"D,W/ S WEDGE	1
F3	ST-ESSSPP2018M	STEALTH,CUSHION,ASSY,ESSENCE,SKIN PROTECTION AND POSITIONING,20"W,18"D,W/ M WEDGE	1
G1	ST-ESSSPP2020	STEALTH,CUSHION,ASSY,ESSENCE,SKIN PROTECTION AND POSITIONING,20"W,20"D	1
G2	ST-ESSSPP2020S	STEALTH,CUSHION,ASSY,ESSENCE,SKIN PROTECTION AND POSITIONING,20"W,20"D,W/ S WEDGE	1
G3	ST-ESSSPP2020M	STEALTH,CUSHION,ASSY,ESSENCE,SKIN PROTECTION AND POSITIONING,20"W,20"D,W/ M WEDGE	1

1. Cushion Assy's

Stealth Simplicity G Complete Cushion Assy	
Size	Ref#
16"W x 16"D	A1
16"W x 18"D	B1
18"W x 16"D	C1
18"W x 18"D	D1
18"W x 20"D	E1
20"W x 16"D	F1
20"W x 18"D	G1
20"W x 20"D	H1
22"W x 18"D	I1
22"W x 20"D	J1
24"W x 18"D	K1





# Stealth Simplicity G Complete Cushion Assy




Item	Part No.	Name	Qty
A1	ST-SIMG1616	STEALTH,CUSHION,ASSY,SIMPLICITY,GENERAL USE,16"W,16"D	1
B1	ST-SIMG1618	STEALTH,CUSHION,ASSY,SIMPLICITY,GENERAL USE,16"W,18"D	1
C1	ST-SIMG1816	STEALTH,CUSHION,ASSY,SIMPLICITY,GENERAL USE,18"W,16"D	1
D1	ST-SIMG1818	STEALTH,CUSHION,ASSY,SIMPLICITY,GENERAL USE,18"W,18"D	1
E1	ST-SIMG1820	STEALTH,CUSHION,ASSY,SIMPLICITY,GENERAL USE,18"W,20"D	1
F1	ST-SIMG2016	STEALTH,CUSHION,ASSY,SIMPLICITY,GENERAL USE,20"W,16"D	1
G1	ST-SIMG2018	STEALTH,CUSHION,ASSY,SIMPLICITY,GENERAL USE,20"W,18"D	1
H1	ST-SIMG2020	STEALTH,CUSHION,ASSY,SIMPLICITY,GENERAL USE,20"W,20"D	1
I1	ST-SIMG2218	STEALTH,CUSHION,ASSY,SIMPLICITY,GENERAL USE,22"W,18"D	1
J1	ST-SIMG2220	STEALTH,CUSHION,ASSY,SIMPLICITY,GENERAL USE,22"W,20"D	1
K1	ST-SIMG2418	STEALTH,CUSHION,ASSY,SIMPLICITY,GENERAL USE,24"W,18"D	1

1. Cushion Assy's

**Stealth Solution Complete Cushion Assy  
with 4-way Stretch Outer Cover**

Size	Ref#	Size	Ref#	Size	Ref#
10"W x 10"D	A1a	16"W x 20"D	A1m	19"W x 19"D	A1y
10"W x 12"D	A1b	17"W x 16"D	A1n	19"W x 20"D	A1z
12"W x 10"D	A1c	17"W x 17"D	A1o	20"W x 16"D	A1aa
12"W x 12"D	A1d	17"W x 18"D	A1p	20"W x 18"D	A1bb
14"W x 12"D	A1e	17"W x 19"D	A1q	20"W x 19"D	A1cc
14"W x 14"D	A1f	17"W x 20"D	A1r	20"W x 20"D	A1dd
14"W x 16"D	A1g	18"W x 16"D	A1s	20"W x 22"D	A1ee
14"W x 18"D	A1h	18"W x 18"D	A1t	22"W x 18"D	A1ff
16"W x 16"D	A1i	18"W x 19"D	A1u	22"W x 20"D	A1gg
16"W x 17"D	A1j	18"W x 20"D	A1v	22"W x 22"D	A1hh
16"W x 18"D	A1k	18"W x 22"D	A1w	24"W x 20"D	A1ii
16"W x 19"D	A1l	19"W x 18"D	A1x	24"W x 22"D	A1jj

\*The Solution cushion assembly comes standard with a urethane inner liner.  
\*\*Optional ABS pan not included and must be selected separately if needed.




Item	Part No.	Name	Qty
A1a	ST-SOLSPP1010	STEALTH,CUSHION,ASSY,SOLUTION,SKIN PROTECTION AND POSITIONING,10"W,10"D	1
A1b	ST-SOLSPP1012	STEALTH,CUSHION,ASSY,SOLUTION,SKIN PROTECTION AND POSITIONING,10"W,12"D	1
A1c	ST-SOLSPP1210	STEALTH,CUSHION,ASSY,SOLUTION,SKIN PROTECTION AND POSITIONING,12"W,10"D	1
A1d	ST-SOLSPP1212	STEALTH,CUSHION,ASSY,SOLUTION,SKIN PROTECTION AND POSITIONING,12"W,12"D	1
A1e	ST-SOLSPP1412	STEALTH,CUSHION,ASSY,SOLUTION,SKIN PROTECTION AND POSITIONING,14"W,12"D	1
A1f	ST-SOLSPP1414	STEALTH,CUSHION,ASSY,SOLUTION,SKIN PROTECTION AND POSITIONING,14"W,14"D	1
A1g	ST-SOLSPP1416	STEALTH,CUSHION,ASSY,SOLUTION,SKIN PROTECTION AND POSITIONING,14"W,16"D	1
A1h	ST-SOLSPP1418	STEALTH,CUSHION,ASSY,SOLUTION,SKIN PROTECTION AND POSITIONING,14"W,18"D	1
A1i	ST-SOLSPP1616	STEALTH,CUSHION,ASSY,SOLUTION,SKIN PROTECTION AND POSITIONING,16"W,16"D	1
A1j	ST-SOLSPP1617	STEALTH,CUSHION,ASSY,SOLUTION,SKIN PROTECTION AND POSITIONING,16"W,17"D	1
A1k	ST-SOLSPP1618	STEALTH,CUSHION,ASSY,SOLUTION,SKIN PROTECTION AND POSITIONING,16"W,18"D	1
A1l	ST-SOLSPP1619	STEALTH,CUSHION,ASSY,SOLUTION,SKIN PROTECTION AND POSITIONING,16"W,19"D	1
A1m	ST-SOLSPP1620	STEALTH,CUSHION,ASSY,SOLUTION,SKIN PROTECTION AND POSITIONING,16"W,20"D	1
A1n	ST-SOLSPP1716	STEALTH,CUSHION,ASSY,SOLUTION,SKIN PROTECTION AND POSITIONING,17"W,16"D	1
A1o	ST-SOLSPP1717	STEALTH,CUSHION,ASSY,SOLUTION,SKIN PROTECTION AND POSITIONING,17"W,17"D	1
A1p	ST-SOLSPP1718	STEALTH,CUSHION,ASSY,SOLUTION,SKIN PROTECTION AND POSITIONING,17"W,18"D	1
A1q	ST-SOLSPP1719	STEALTH,CUSHION,ASSY,SOLUTION,SKIN PROTECTION AND POSITIONING,17"W,19"D	1
A1r	ST-SOLSPP1720	STEALTH,CUSHION,ASSY,SOLUTION,SKIN PROTECTION AND POSITIONING,17"W,20"D	1
A1s	ST-SOLSPP1816	STEALTH,CUSHION,ASSY,SOLUTION,SKIN PROTECTION AND POSITIONING,18"W,16"D	1
A1t	ST-SOLSPP1818	STEALTH,CUSHION,ASSY,SOLUTION,SKIN PROTECTION AND POSITIONING,18"W,18"D	1
A1u	ST-SOLSPP1819	STEALTH,CUSHION,ASSY,SOLUTION,SKIN PROTECTION AND POSITIONING,18"W,19"D	1
A1v	ST-SOLSPP1820	STEALTH,CUSHION,ASSY,SOLUTION,SKIN PROTECTION AND POSITIONING,18"W,20"D	1
A1w	ST-SOLSPP1822	STEALTH,CUSHION,ASSY,SOLUTION,SKIN PROTECTION AND POSITIONING,18"W,22"D	1
A1x	ST-SOLSPP1918	STEALTH,CUSHION,ASSY,SOLUTION,SKIN PROTECTION AND POSITIONING,19"W,18"D	1
A1y	ST-SOLSPP1919	STEALTH,CUSHION,ASSY,SOLUTION,SKIN PROTECTION AND POSITIONING,19"W,19"D	1
A1z	ST-SOLSPP1920	STEALTH,CUSHION,ASSY,SOLUTION,SKIN PROTECTION AND POSITIONING,19"W,20"D	1
A1a	ST-SOLSPP2016	STEALTH,CUSHION,ASSY,SOLUTION,SKIN PROTECTION AND POSITIONING,20"W,16"D	1
A1b	ST-SOLSPP2018	STEALTH,CUSHION,ASSY,SOLUTION,SKIN PROTECTION AND POSITIONING,20"W,18"D	1
A1c	ST-SOLSPP2019	STEALTH,CUSHION,ASSY,SOLUTION,SKIN PROTECTION AND POSITIONING,20"W,19"D	1
A1d	ST-SOLSPP2020	STEALTH,CUSHION,ASSY,SOLUTION,SKIN PROTECTION AND POSITIONING,20"W,20"D	1
A1e	ST-SOLSPP2022	STEALTH,CUSHION,ASSY,SOLUTION,SKIN PROTECTION AND POSITIONING,20"W,22"D	1
A1ff	ST-SOLSPP2218	STEALTH,CUSHION,ASSY,SOLUTION,SKIN PROTECTION AND POSITIONING,22"W,18"D	1

### Stealth Solution Complete Cushion Assy with 4-way Stretch Outer Cover

Size	Ref#	Size	Ref#	Size	Ref#
10"W x 10"D	A1a	16"W x 20"D	A1m	19"W x 19"D	A1y
10"W x 12"D	A1b	17"W x 16"D	A1n	19"W x 20"D	A1z
12"W x 10"D	A1c	17"W x 17"D	A1o	20"W x 16"D	A1aa
12"W x 12"D	A1d	17"W x 18"D	A1p	20"W x 18"D	A1bb
14"W x 12"D	A1e	17"W x 19"D	A1q	20"W x 19"D	A1cc
14"W x 14"D	A1f	17"W x 20"D	A1r	20"W x 20"D	A1dd
14"W x 16"D	A1g	18"W x 16"D	A1s	20"W x 22"D	A1ee
14"W x 18"D	A1h	18"W x 18"D	A1t	22"W x 18"D	A1ff
16"W x 16"D	A1i	18"W x 19"D	A1u	22"W x 20"D	A1gg
16"W x 17"D	A1j	18"W x 20"D	A1v	22"W x 22"D	A1hh
16"W x 18"D	A1k	18"W x 22"D	A1w	24"W x 20"D	A1ii
16"W x 19"D	A1l	19"W x 18"D	A1x	24"W x 22"D	A1jj

\*The Solution cushion assembly comes standard with a urethane inner liner.  
 \*\*Optional ABS pan not included and must be selected separately if needed.




Item	Part No.	Name	Qty
A1g	ST-SOLSPP2220	STEALTH,CUSHION,ASSY,SOLUTION,SKIN PROTECTION AND POSITIONING,22"W,20"D	1
A1h	ST-SOLSPP2222	STEALTH,CUSHION,ASSY,SOLUTION,SKIN PROTECTION AND POSITIONING,22"W,22"D	1
A1ii	ST-SOLSPP2420	STEALTH,CUSHION,ASSY,SOLUTION,SKIN PROTECTION AND POSITIONING,24"W,20"D	1
A1jj	ST-SOLSPP2422	STEALTH,CUSHION,ASSY,SOLUTION,SKIN PROTECTION AND POSITIONING,24"W,22"D	1

1. Cushion Assy's

Stealth Spectrum Foam Complete Cushion Assy with 4-way Stretch Outer Cover							
Size	Ref#	Size	Ref#	Size	Ref#	Size	Ref#
16"W x 16"D	A1a	20"W x 16"D	A1m	24"W x 22"D	A1y	30"W x 20"D	A1kk
16"W x 18"D	A1b	20"W x 18"D	A1n	24"W x 24"D	A1z	30"W x 22"D	A1ll
16"W x 20"D	A1c	20"W x 20"D	A1o	26"W x 18"D	A1aa	30"W x 24"D	A1mm
17"W x 17"D	A1d	20"W x 22"D	A1p	26"W x 20"D	A1bb	30"W x 28"D	A1nn
17"W x 18"D	A1e	20"W x 24"D	A1q	26"W x 22"D	A1cc	32"W x 26"D	A1oo
17"W x 19"D	A1f	22"W x 18"D	A1r	26"W x 24"D	A1dd	32"W x 28"D	A1pp
18"W x 16"D	A1g	22"W x 20"D	A1s	28"W x 18"D	A1ee		
18"W x 18"D	A1h	22"W x 22"D	A1t	28"W x 20"D	A1ff		
18"W x 20"D	A1i	22"W x 24"D	A1u	28"W x 22"D	A1gg		
18"W x 22"D	A1j	24"W x 16"D	A1v	28"W x 24"D	A1hh		
19"W x 17"D	A1k	24"W x 18"D	A1w	28"W x 28"D	A1ii		
19"W x 19"D	A1l	24"W x 20"D	A1x	30"W x 18"D	A1jj		

\*The Spectrum Foam cushion assembly comes standard with a urethane inner liner.  
 \*\*Optional ABS pan not included and must be selected separately if needed.




Item	Part No.	Name	Qty
A1a	ST-SPCFSP1616	STEALTH,CUSHION,ASSY,SPECTRUM FOAM,SKIN PROTECTION AND POSITIONING,16"W,16"D	1
A1b	ST-SPCFSP1618	STEALTH,CUSHION,ASSY,SPECTRUM FOAM,SKIN PROTECTION AND POSITIONING,16"W,18"D	1
A1c	ST-SPCFSP1620	STEALTH,CUSHION,ASSY,SPECTRUM FOAM,SKIN PROTECTION AND POSITIONING,16"W,20"D	1
A1d	ST-SPCFSP1717	STEALTH,CUSHION,ASSY,SPECTRUM FOAM,SKIN PROTECTION AND POSITIONING,17"W,17"D	1
A1e	ST-SPCFSP1718	STEALTH,CUSHION,ASSY,SPECTRUM FOAM,SKIN PROTECTION AND POSITIONING,17"W,18"D	1
A1f	ST-SPCFSP1719	STEALTH,CUSHION,ASSY,SPECTRUM FOAM,SKIN PROTECTION AND POSITIONING,17"W,19"D	1
A1g	ST-SPCFSP1816	STEALTH,CUSHION,ASSY,SPECTRUM FOAM,SKIN PROTECTION AND POSITIONING,18"W,16"D	1
A1h	ST-SPCFSP1818	STEALTH,CUSHION,ASSY,SPECTRUM FOAM,SKIN PROTECTION AND POSITIONING,18"W,18"D	1
A1i	ST-SPCFSP1820	STEALTH,CUSHION,ASSY,SPECTRUM FOAM,SKIN PROTECTION AND POSITIONING,18"W,20"D	1
A1j	ST-SPCFSP1822	STEALTH,CUSHION,ASSY,SPECTRUM FOAM,SKIN PROTECTION AND POSITIONING,18"W,22"D	1
A1k	ST-SPCFSP1917	STEALTH,CUSHION,ASSY,SPECTRUM FOAM,SKIN PROTECTION AND POSITIONING,19"W,17"D	1
A1l	ST-SPCFSP1919	STEALTH,CUSHION,ASSY,SPECTRUM FOAM,SKIN PROTECTION AND POSITIONING,19"W,19"D	1
A1m	ST-SPCFSP2016	STEALTH,CUSHION,ASSY,SPECTRUM FOAM,SKIN PROTECTION AND POSITIONING,20"W,16"D	1
A1n	ST-SPCFSP2018	STEALTH,CUSHION,ASSY,SPECTRUM FOAM,SKIN PROTECTION AND POSITIONING,20"W,18"D	1
A1o	ST-SPCFSP2020	STEALTH,CUSHION,ASSY,SPECTRUM FOAM,SKIN PROTECTION AND POSITIONING,20"W,20"D	1
A1p	ST-SPCFSP2022	STEALTH,CUSHION,ASSY,SPECTRUM FOAM,SKIN PROTECTION AND POSITIONING,20"W,22"D	1
A1q	ST-SPCFSP2024	STEALTH,CUSHION,ASSY,SPECTRUM FOAM,SKIN PROTECTION AND POSITIONING,20"W,24"D	1
A1r	ST-SPCFSP2218	STEALTH,CUSHION,ASSY,SPECTRUM FOAM,SKIN PROTECTION AND POSITIONING,22"W,18"D	1
A1s	ST-SPCFSP2220	STEALTH,CUSHION,ASSY,SPECTRUM FOAM,SKIN PROTECTION AND POSITIONING,22"W,20"D	1
A1t	ST-SPCFSP2222	STEALTH,CUSHION,ASSY,SPECTRUM FOAM,SKIN PROTECTION AND POSITIONING,22"W,22"D	1
A1u	ST-SPCFSP2224	STEALTH,CUSHION,ASSY,SPECTRUM FOAM,SKIN PROTECTION AND POSITIONING,22"W,24"D	1
A1v	ST-SPCFSP2416	STEALTH,CUSHION,ASSY,SPECTRUM FOAM,SKIN PROTECTION AND POSITIONING,24"W,16"D	1
A1w	ST-SPCFSP2418	STEALTH,CUSHION,ASSY,SPECTRUM FOAM,SKIN PROTECTION AND POSITIONING,24"W,18"D	1
A1x	ST-SPCFSP2420	STEALTH,CUSHION,ASSY,SPECTRUM FOAM,SKIN PROTECTION AND POSITIONING,24"W,20"D	1
A1y	ST-SPCFSP2422	STEALTH,CUSHION,ASSY,SPECTRUM FOAM,SKIN PROTECTION AND POSITIONING,24"W,22"D	1
A1z	ST-SPCFSP2424	STEALTH,CUSHION,ASSY,SPECTRUM FOAM,SKIN PROTECTION AND POSITIONING,24"W,24"D	1
A1a	ST-SPCFSP2618	STEALTH,CUSHION,ASSY,SPECTRUM FOAM,SKIN PROTECTION AND POSITIONING,26"W,18"D	1
A1b	ST-SPCFSP2620	STEALTH,CUSHION,ASSY,SPECTRUM FOAM,SKIN PROTECTION AND POSITIONING,26"W,20"D	1
A1c	ST-SPCFSP2622	STEALTH,CUSHION,ASSY,SPECTRUM FOAM,SKIN PROTECTION AND POSITIONING,26"W,22"D	1
A1d	ST-SPCFSP2624	STEALTH,CUSHION,ASSY,SPECTRUM FOAM,SKIN PROTECTION AND POSITIONING,26"W,24"D	1
A1e	ST-SPCFSP2818	STEALTH,CUSHION,ASSY,SPECTRUM FOAM,SKIN PROTECTION AND POSITIONING,28"W,18"D	1
A1ff	ST-SPCFSP2820	STEALTH,CUSHION,ASSY,SPECTRUM FOAM,SKIN PROTECTION AND POSITIONING,28"W,20"D	1

### Stealth Spectrum Foam Complete Cushion Assy with 4-way Stretch Outer Cover

Size	Ref#	Size	Ref#	Size	Ref#	Size	Ref#
16"W x 16"D	A1a	20"W x 16"D	A1m	24"W x 22"D	A1y	30"W x 20"D	A1kk
16"W x 18"D	A1b	20"W x 18"D	A1n	24"W x 24"D	A1z	30"W x 22"D	A1ll
16"W x 20"D	A1c	20"W x 20"D	A1o	26"W x 18"D	A1aa	30"W x 24"D	A1mm
17"W x 17"D	A1d	20"W x 22"D	A1p	26"W x 20"D	A1bb	30"W x 28"D	A1nn
17"W x 18"D	A1e	20"W x 24"D	A1q	26"W x 22"D	A1cc	32"W x 26"D	A1oo
17"W x 19"D	A1f	22"W x 18"D	A1r	26"W x 24"D	A1dd	32"W x 28"D	A1pp
18"W x 16"D	A1g	22"W x 20"D	A1s	28"W x 18"D	A1ee		
18"W x 18"D	A1h	22"W x 22"D	A1t	28"W x 20"D	A1ff		
18"W x 20"D	A1i	22"W x 24"D	A1u	28"W x 22"D	A1gg		
18"W x 22"D	A1j	24"W x 16"D	A1v	28"W x 24"D	A1hh		
19"W x 17"D	A1k	24"W x 18"D	A1w	28"W x 28"D	A1ii		
19"W x 19"D	A1l	24"W x 20"D	A1x	30"W x 18"D	A1jj		

\*The Spectrum Foam cushion assembly comes standard with a urethane inner liner.  
 \*\*Optional ABS pan not included and must be selected separately if needed.






Item	Part No.	Name	Qty
A1g	ST-SPCFSP2822	STEALTH,CUSHION,ASSY,SPECTRUM FOAM,SKIN PROTECTION AND POSITIONING,28"W,22"D	1
A1h	ST-SPCFSP2824	STEALTH,CUSHION,ASSY,SPECTRUM FOAM,SKIN PROTECTION AND POSITIONING,28"W,24"D	1
A1ii	ST-SPCFSP2828	STEALTH,CUSHION,ASSY,SPECTRUM FOAM,SKIN PROTECTION AND POSITIONING,28"W,28"D	1
A1jj	ST-SPCFSP3018	STEALTH,CUSHION,ASSY,SPECTRUM FOAM,SKIN PROTECTION AND POSITIONING,30"W,18"D	1
A1kk	ST-SPCFSP3020	STEALTH,CUSHION,ASSY,SPECTRUM FOAM,SKIN PROTECTION AND POSITIONING,30"W,20"D	1
A1ll	ST-SPCFSP3022	STEALTH,CUSHION,ASSY,SPECTRUM FOAM,SKIN PROTECTION AND POSITIONING,30"W,22"D	1
A1mm	ST-SPCFSP3024	STEALTH,CUSHION,ASSY,SPECTRUM FOAM,SKIN PROTECTION AND POSITIONING,30"W,24"D	1
A1nn	ST-SPCFSP3028	STEALTH,CUSHION,ASSY,SPECTRUM FOAM,SKIN PROTECTION AND POSITIONING,30"W,28"D	1
A1oo	ST-SPCFSP3226	STEALTH,CUSHION,ASSY,SPECTRUM FOAM,SKIN PROTECTION AND POSITIONING,32"W,26"D	1
A1pp	ST-SPCFSP3228	STEALTH,CUSHION,ASSY,SPECTRUM FOAM,SKIN PROTECTION AND POSITIONING,32"W,28"D	1

1. Cushion Assy's

Stealth Spectrum Gel Complete Cushion Assy					
Size	Ref#	Size	Ref#	Size	Ref#
14"W x 18"D	A1a	18"W x 20"D	A1n	22"W x 16"D	A1aa
16"W x 16"D	A1b	18"W x 22"D	A1o	22"W x 18"D	A1bb
16"W x 18"D	A1c	18"W x 24"D	A1p	22"W x 19"D	A1cc
16"W x 19"D	A1d	19"W x 17"D	A1q	22"W x 20"D	A1dd
16"W x 20"D	A1e	19"W x 18"D	A1r	22"W x 22"D	A1ee
17"W x 16"D	A1f	19"W x 19"D	A1s	22"W x 24"D	A1ff
17"W x 17"D	A1g	19"W x 20"D	A1t	24"W x 16"D	A1gg
17"W x 18"D	A1h	20"W x 16"D	A1u	24"W x 18"D	A1hh
17"W x 19"D	A1i	20"W x 18"D	A1v	24"W x 20"D	A1ii
17"W x 20"D	A1j	20"W x 19"D	A1w	24"W x 22"D	A1jj
18"W x 16"D	A1k	20"W x 20"D	A1x	24"W x 24"D	A1kk
18"W x 17"D	A1l	20"W x 22"D	A1y		
18"W x 18"D	A1m	20"W x 24"D	A1z		

\*The Spectrum Gel cushion assembly comes standard with a urethane inner liner.  
 \*\*Optional ABS pan not included and must be selected separately if needed.




Item	Part No.	Name	Qty
A1a	ST-SPCGSPP1418	STEALTH,CUSHION,ASSY,SPECTRUM GEL,SKIN PROTECTION AND POSITIONING,14"W,18"D	1
A1b	ST-SPCGSPP1616	STEALTH,CUSHION,ASSY,SPECTRUM GEL,SKIN PROTECTION AND POSITIONING,16"W,16"D	1
A1c	ST-SPCGSPP1618	STEALTH,CUSHION,ASSY,SPECTRUM GEL,SKIN PROTECTION AND POSITIONING,16"W,18"D	1
A1d	ST-SPCGSPP1619	STEALTH,CUSHION,ASSY,SPECTRUM GEL,SKIN PROTECTION AND POSITIONING,16"W,19"D	1
A1e	ST-SPCGSPP1620	STEALTH,CUSHION,ASSY,SPECTRUM GEL,SKIN PROTECTION AND POSITIONING,16"W,20"D	1
A1f	ST-SPCGSPP1716	STEALTH,CUSHION,ASSY,SPECTRUM GEL,SKIN PROTECTION AND POSITIONING,17"W,16"D	1
A1g	ST-SPCGSPP1717	STEALTH,CUSHION,ASSY,SPECTRUM GEL,SKIN PROTECTION AND POSITIONING,17"W,17"D	1
A1h	ST-SPCGSPP1718	STEALTH,CUSHION,ASSY,SPECTRUM GEL,SKIN PROTECTION AND POSITIONING,17"W,18"D	1
A1i	ST-SPCGSPP1719	STEALTH,CUSHION,ASSY,SPECTRUM GEL,SKIN PROTECTION AND POSITIONING,17"W,19"D	1
A1j	ST-SPCGSPP1720	STEALTH,CUSHION,ASSY,SPECTRUM GEL,SKIN PROTECTION AND POSITIONING,17"W,20"D	1
A1k	ST-SPCGSPP1816	STEALTH,CUSHION,ASSY,SPECTRUM GEL,SKIN PROTECTION AND POSITIONING,18"W,16"D	1
A1l	ST-SPCGSPP1817	STEALTH,CUSHION,ASSY,SPECTRUM GEL,SKIN PROTECTION AND POSITIONING,18"W,17"D	1
A1m	ST-SPCGSPP1818	STEALTH,CUSHION,ASSY,SPECTRUM GEL,SKIN PROTECTION AND POSITIONING,18"W,18"D	1
A1n	ST-SPCGSPP1820	STEALTH,CUSHION,ASSY,SPECTRUM GEL,SKIN PROTECTION AND POSITIONING,18"W,20"D	1
A1o	ST-SPCGSPP1822	STEALTH,CUSHION,ASSY,SPECTRUM GEL,SKIN PROTECTION AND POSITIONING,18"W,22"D	1
A1p	ST-SPCGSPP1824	STEALTH,CUSHION,ASSY,SPECTRUM GEL,SKIN PROTECTION AND POSITIONING,18"W,24"D	1
A1q	ST-SPCGSPP1917	STEALTH,CUSHION,ASSY,SPECTRUM GEL,SKIN PROTECTION AND POSITIONING,19"W,17"D	1
A1r	ST-SPCGSPP1918	STEALTH,CUSHION,ASSY,SPECTRUM GEL,SKIN PROTECTION AND POSITIONING,19"W,18"D	1
A1s	ST-SPCGSPP1919	STEALTH,CUSHION,ASSY,SPECTRUM GEL,SKIN PROTECTION AND POSITIONING,19"W,19"D	1
A1t	ST-SPCGSPP1920	STEALTH,CUSHION,ASSY,SPECTRUM GEL,SKIN PROTECTION AND POSITIONING,19"W,20"D	1
A1u	ST-SPCGSPP2016	STEALTH,CUSHION,ASSY,SPECTRUM GEL,SKIN PROTECTION AND POSITIONING,20"W,16"D	1
A1v	ST-SPCGSPP2018	STEALTH,CUSHION,ASSY,SPECTRUM GEL,SKIN PROTECTION AND POSITIONING,20"W,18"D	1
A1w	ST-SPCGSPP2019	STEALTH,CUSHION,ASSY,SPECTRUM GEL,SKIN PROTECTION AND POSITIONING,20"W,19"D	1
A1x	ST-SPCGSPP2020	STEALTH,CUSHION,ASSY,SPECTRUM GEL,SKIN PROTECTION AND POSITIONING,20"W,20"D	1
A1y	ST-SPCGSPP2022	STEALTH,CUSHION,ASSY,SPECTRUM GEL,SKIN PROTECTION AND POSITIONING,20"W,22"D	1
A1z	ST-SPCGSPP2024	STEALTH,CUSHION,ASSY,SPECTRUM GEL,SKIN PROTECTION AND POSITIONING,20"W,24"D	1
A1a a	ST-SPCGSPP2216	STEALTH,CUSHION,ASSY,SPECTRUM GEL,SKIN PROTECTION AND POSITIONING,22"W,16"D	1
A1b b	ST-SPCGSPP2218	STEALTH,CUSHION,ASSY,SPECTRUM GEL,SKIN PROTECTION AND POSITIONING,22"W,18"D	1
A1c c	ST-SPCGSPP2219	STEALTH,CUSHION,ASSY,SPECTRUM GEL,SKIN PROTECTION AND POSITIONING,22"W,19"D	1
A1d d	ST-SPCGSPP2220	STEALTH,CUSHION,ASSY,SPECTRUM GEL,SKIN PROTECTION AND POSITIONING,22"W,20"D	1
A1e e	ST-SPCGSPP2222	STEALTH,CUSHION,ASSY,SPECTRUM GEL,SKIN PROTECTION AND POSITIONING,22"W,22"D	1
A1ff	ST-SPCGSPP2224	STEALTH,CUSHION,ASSY,SPECTRUM GEL,SKIN PROTECTION AND POSITIONING,22"W,24"D	1

### Stealth Spectrum Gel Complete Cushion Assy

Size	Ref#	Size	Ref#	Size	Ref#
14"W x 18"D	A1a	18"W x 20"D	A1n	22"W x 16"D	A1aa
16"W x 16"D	A1b	18"W x 22"D	A1o	22"W x 18"D	A1bb
16"W x 18"D	A1c	18"W x 24"D	A1p	22"W x 19"D	A1cc
16"W x 19"D	A1d	19"W x 17"D	A1q	22"W x 20"D	A1dd
16"W x 20"D	A1e	19"W x 18"D	A1r	22"W x 22"D	A1ee
17"W x 16"D	A1f	19"W x 19"D	A1s	22"W x 24"D	A1ff
17"W x 17"D	A1g	19"W x 20"D	A1t	24"W x 16"D	A1gg
17"W x 18"D	A1h	20"W x 16"D	A1u	24"W x 18"D	A1hh
17"W x 19"D	A1i	20"W x 18"D	A1v	24"W x 20"D	A1ii
17"W x 20"D	A1j	20"W x 19"D	A1w	24"W x 22"D	A1jj
18"W x 16"D	A1k	20"W x 20"D	A1x	24"W x 24"D	A1kk
18"W x 17"D	A1l	20"W x 22"D	A1y		
18"W x 18"D	A1m	20"W x 24"D	A1z		

\*The Spectrum Gel cushion assembly comes standard with a urethane inner liner.

\*\*Optional ABS pan not included and must be selected separately if needed.



Item	Part No.	Name	Qty
A1g	ST-SPCGSPP2416	STEALTH,CUSHION,ASSY,SPECTRUM GEL,SKIN PROTECTION AND POSITIONING,24"W,16"D	1
A1h	ST-SPCGSPP2418	STEALTH,CUSHION,ASSY,SPECTRUM GEL,SKIN PROTECTION AND POSITIONING,24"W,18"D	1
A1ii	ST-SPCGSPP2420	STEALTH,CUSHION,ASSY,SPECTRUM GEL,SKIN PROTECTION AND POSITIONING,24"W,20"D	1
A1jj	ST-SPCGSPP2422	STEALTH,CUSHION,ASSY,SPECTRUM GEL,SKIN PROTECTION AND POSITIONING,24"W,22"D	1
A1kk	ST-SPCGSPP2424	STEALTH,CUSHION,ASSY,SPECTRUM GEL,SKIN PROTECTION AND POSITIONING,24"W,24"D	1

## 1. Cushion Assy's

### Stealth TRU-Comfort 2 Complete Cushion Assy

Size	Ref#	Size	Ref#
16"W x 14"D	A1a	20"W x 18"D	A1j
16"W x 16"D	A1b	20"W x 20"D	A1k
16"W x 18"D	A1c	20"W x 22"D	A1l
16"W x 20"D	A1d	20"W x 24"D	A1m
18"W x 16"D	A1e	22"W x 18"D	A1n
18"W x 18"D	A1f	22"W x 20"D	A1o
18"W x 20"D	A1g	22"W x 22"D	A1p
18"W x 22"D	A1h	22"W x 24"D	A1q
20"W x 16"D	A1i		



Item	Part No.	Name	Qty
A1a	ST-TC2SPP1614	STEALTH,CUSHION,ASSY,TRU-COMFORT 2,SKIN PROTECTION AND POSITIONING,16"W,14"D	1
A1b	ST-TC2SPP1616	STEALTH,CUSHION,ASSY,TRU-COMFORT 2,SKIN PROTECTION AND POSITIONING,16"W,16"D	1
A1c	ST-TC2SPP1618	STEALTH,CUSHION,ASSY,TRU-COMFORT 2,SKIN PROTECTION AND POSITIONING,16"W,18"D	1
A1d	ST-TC2SPP1620	STEALTH,CUSHION,ASSY,TRU-COMFORT 2,SKIN PROTECTION AND POSITIONING,16"W,20"D	1
A1e	ST-TC2SPP1816	STEALTH,CUSHION,ASSY,TRU-COMFORT 2,SKIN PROTECTION AND POSITIONING,18"W,16"D	1
A1f	ST-TC2SPP1818	STEALTH,CUSHION,ASSY,TRU-COMFORT 2,SKIN PROTECTION AND POSITIONING,18"W,18"D	1
A1g	ST-TC2SPP1820	STEALTH,CUSHION,ASSY,TRU-COMFORT 2,SKIN PROTECTION AND POSITIONING,18"W,20"D	1
A1h	ST-TC2SPP1822	STEALTH,CUSHION,ASSY,TRU-COMFORT 2,SKIN PROTECTION AND POSITIONING,18"W,22"D	1
A1i	ST-TC2SPP2016	STEALTH,CUSHION,ASSY,TRU-COMFORT 2,SKIN PROTECTION AND POSITIONING,20"W,16"D	1
A1j	ST-TC2SPP2018	STEALTH,CUSHION,ASSY,TRU-COMFORT 2,SKIN PROTECTION AND POSITIONING,20"W,18"D	1
A1k	ST-TC2SPP2020	STEALTH,CUSHION,ASSY,TRU-COMFORT 2,SKIN PROTECTION AND POSITIONING,20"W,20"D	1
A1l	ST-TC2SPP2022	STEALTH,CUSHION,ASSY,TRU-COMFORT 2,SKIN PROTECTION AND POSITIONING,20"W,22"D	1
A1m	ST-TC2SPP2024	STEALTH,CUSHION,ASSY,TRU-COMFORT 2,SKIN PROTECTION AND POSITIONING,20"W,24"D	1
A1n	ST-TC2SPP2218	STEALTH,CUSHION,ASSY,TRU-COMFORT 2,SKIN PROTECTION AND POSITIONING,22"W,18"D	1
A1o	ST-TC2SPP2220	STEALTH,CUSHION,ASSY,TRU-COMFORT 2,SKIN PROTECTION AND POSITIONING,22"W,20"D	1
A1p	ST-TC2SPP2222	STEALTH,CUSHION,ASSY,TRU-COMFORT 2,SKIN PROTECTION AND POSITIONING,22"W,22"D	1
A1q	ST-TC2SPP2224	STEALTH,CUSHION,ASSY,TRU-COMFORT 2,SKIN PROTECTION AND POSITIONING,22"W,24"D	1

2. Outer Covers




\*Please note: Shown without wedges



Item	Part No.	Name	Qty
A1	ESSSPP-NBLK1616	CUSHION COVER,ASSY,ESS,4 WAY STRETCH,16"W,16"D,W/PVC	1
B1	ESSSPP-NBLK1618	CUSHION COVER,ASSY,ESS,4 WAY STRETCH,16"W,18"D,W/PVC	1
C1	ESSSPP-NBLK1816	CUSHION COVER,ASSY,ESS,4 WAY STRETCH,18"W,16"D,W/PVC	1
D1	ESSSPP-NBLK1818	CUSHION COVER,ASSY,ESS,4 WAY STRETCH,18"W,18"D,W/PVC	1
E1	ESSSPP-NBLK1820	CUSHION COVER,ASSY,ESS,4 WAY STRETCH,18"W,20"D,W/PVC	1
F1	ESSSPP-NBLK2018	CUSHION COVER,ASSY,ESS,4 WAY STRETCH,20"W,18"D,W/PVC	1
G1	ESSSPP-NBLK2020	CUSHION COVER,ASSY,ESS,4 WAY STRETCH,20"W,20"D,W/PVC	1


2. Outer Covers



\*Please note: Shown without wedges

Item	Part No.	Name	Qty
A1	ESSPLSPP-NBLK1616	CUSHION COVER,ASSY,ESSPL,4 WAY STRETCH,16"W,16"D,W/PVC	1
B1	ESSPLSPP-NBLK1618	CUSHION COVER,ASSY,ESSPL,4 WAY STRETCH,16"W,18"D,W/PVC	1
C1	ESSPLSPP-NBLK1816	CUSHION COVER,ASSY,ESSPL,4 WAY STRETCH,18"W,16"D,W/PVC	1
D1	ESSPLSPP-NBLK1818	CUSHION COVER,ASSY,ESSPL,4 WAY STRETCH,18"W,18"D,W/PVC	1
E1	ESSPLSPP-NBLK1820	CUSHION COVER,ASSY,ESSPL,4 WAY STRETCH,18"W,20"D,W/PVC	1
F1	ESSPLSPP-NBLK2018	CUSHION COVER,ASSY,ESSPL,4 WAY STRETCH,20"W,18"D,W/PVC	1
G1	ESSPLSPP-NBLK2020	CUSHION COVER,ASSY,ESSPL,4 WAY STRETCH,20"W,20"D,W/PVC	1

2. Outer Covers




Stealth Simplicity G Tek2 Outer Cover	
Size	Ref#
16"W x 16"D	A1
16"W x 18"D	B1
18"W x 16"D	C1
18"W x 18"D	D1
18"W x 20"D	E1
20"W x 16"D	F1
20"W x 18"D	G1
20"W x 20"D	H1
22"W x 18"D	I1
22"W x 20"D	J1
24"W x 18"D	K1

## Stealth Simplicity G Tek2 Outer Cushion Cover

Item	Part No.	Name	Qty
A1	SIMG-TBLK1616	CUSHION COVER,ASSY,SIM,2 WAY STRETCH,16"W,16"D,W/PVC	1
B1	SIMG-TBLK1618	CUSHION COVER,ASSY,SIM,2 WAY STRETCH,16"W,18"D,W/PVC	1
C1	SIMG-TBLK1816	CUSHION COVER,ASSY,SIM,2 WAY STRETCH,18"W,16"D,W/PVC	1
D1	SIMG-TBLK1818	CUSHION COVER,ASSY,SIM,2 WAY STRETCH,18"W,18"D,W/PVC	1
E1	SIMG-TBLK1820	CUSHION COVER,ASSY,SIM,2 WAY STRETCH,18"W,20"D,W/PVC	1
F1	SIMG-TBLK2016	CUSHION COVER,ASSY,SIM,2 WAY STRETCH,20"W,16"D,W/PVC	1
G1	SIMG-TBLK2018	CUSHION COVER,ASSY,SIM,2 WAY STRETCH,20"W,18"D,W/PVC	1
H1	SIMG-TBLK2020	CUSHION COVER,ASSY,SIM,2 WAY STRETCH,20"W,20"D,W/PVC	1
I1	SIMG-TBLK2218	CUSHION COVER,ASSY,SIM,2 WAY STRETCH,22"W,18"D,W/PVC	1
J1	SIMG-TBLK2220	CUSHION COVER,ASSY,SIM,2 WAY STRETCH,22"W,20"D,W/PVC	1
K1	SIMG-TBLK2418	CUSHION COVER,ASSY,SIM,2 WAY STRETCH,24"W,18"D,W/PVC	1

2. Outer Covers

Stealth Solution 4-way Stretch Outer Cushion Cover					
Size	Ref#	Size	Ref#	Size	Ref#
10"W x 10"D	A1a	16"W x 20"D	A1m	19"W x 19"D	A1y
10"W x 12"D	A1b	17"W x 16"D	A1n	19"W x 20"D	A1z
12"W x 10"D	A1c	17"W x 17"D	A1o	20"W x 16"D	A1aa
12"W x 12"D	A1d	17"W x 18"D	A1p	20"W x 18"D	A1bb
14"W x 12"D	A1e	17"W x 19"D	A1q	20"W x 19"D	A1cc
14"W x 14"D	A1f	17"W x 20"D	A1r	20"W x 20"D	A1dd
14"W x 16"D	A1g	18"W x 16"D	A1s	20"W x 22"D	A1ee
14"W x 18"D	A1h	18"W x 18"D	A1t	22"W x 18"D	A1ff
16"W x 16"D	A1i	18"W x 19"D	A1u	22"W x 20"D	A1gg
16"W x 17"D	A1j	18"W x 20"D	A1v	22"W x 22"D	A1hh
16"W x 18"D	A1k	18"W x 22"D	A1w	24"W x 20"D	A1ii
16"W x 19"D	A1l	19"W x 18"D	A1x	24"W x 22"D	A1jj




# Stealth Solution 4-way Stretch Outer Cushion Cover



Ite m	Part No.	Name	Qty
A1a	SOLSPP-NBLK1010	CUSHION COVER,ASSY,SOL,4 WAY STRETCH,10"W,10"D,W/PVC	1
A1b	SOLSPP-NBLK1012	CUSHION COVER,ASSY,SOL,4 WAY STRETCH,10"W,12"D,W/PVC	1
A1c	SOLSPP-NBLK1210	CUSHION COVER,ASSY,SOL,4 WAY STRETCH,12"W,10"D,W/PVC	1
A1d	SOLSPP-NBLK1212	CUSHION COVER,ASSY,SOL,4 WAY STRETCH,12"W,12"D,W/PVC	1
A1e	SOLSPP-NBLK1412	CUSHION COVER,ASSY,SOL,4 WAY STRETCH,14"W,12"D,W/PVC	1
A1f	SOLSPP-NBLK1414	CUSHION COVER,ASSY,SOL,4 WAY STRETCH,14"W,14"D,W/PVC	1
A1g	SOLSPP-NBLK1416	CUSHION COVER,ASSY,SOL,4 WAY STRETCH,14"W,16"D,W/PVC	1
A1h	SOLSPP-NBLK1418	CUSHION COVER,ASSY,SOL,4 WAY STRETCH,14"W,18"D,W/PVC	1
A1i	SOLSPP-NBLK1616	CUSHION COVER,ASSY,SOL,4 WAY STRETCH,16"W,16"D,W/PVC	1
A1j	SOLSPP-NBLK1617	CUSHION COVER,ASSY,SOL,4 WAY STRETCH,16"W,17"D,W/PVC	1
A1k	SOLSPP-NBLK1618	CUSHION COVER,ASSY,SOL,4 WAY STRETCH,16"W,18"D,W/PVC	1
A1l	SOLSPP-NBLK1619	CUSHION COVER,ASSY,SOL,4 WAY STRETCH,16"W,19"D,W/PVC	1
A1m	SOLSPP-NBLK1620	CUSHION COVER,ASSY,SOL,4 WAY STRETCH,16"W,20"D,W/PVC	1
A1n	SOLSPP-NBLK1716	CUSHION COVER,ASSY,SOL,4 WAY STRETCH,17"W,16"D,W/PVC	1
A1o	SOLSPP-NBLK1717	CUSHION COVER,ASSY,SOL,4 WAY STRETCH,17"W,17"D,W/PVC	1
A1p	SOLSPP-NBLK1718	CUSHION COVER,ASSY,SOL,4 WAY STRETCH,17"W,18"D,W/PVC	1
A1q	SOLSPP-NBLK1719	CUSHION COVER,ASSY,SOL,4 WAY STRETCH,17"W,19"D,W/PVC	1
A1r	SOLSPP-NBLK1720	CUSHION COVER,ASSY,SOL,4 WAY STRETCH,17"W,20"D,W/PVC	1
A1s	SOLSPP-NBLK1816	CUSHION COVER,ASSY,SOL,4 WAY STRETCH,18"W,16"D,W/PVC	1
A1t	SOLSPP-NBLK1818	CUSHION COVER,ASSY,SOL,4 WAY STRETCH,18"W,18"D,W/PVC	1
A1u	SOLSPP-NBLK1819	CUSHION COVER,ASSY,SOL,4 WAY STRETCH,18"W,19"D,W/PVC	1
A1v	SOLSPP-NBLK1820	CUSHION COVER,ASSY,SOL,4 WAY STRETCH,18"W,20"D,W/PVC	1
A1w	SOLSPP-NBLK1822	CUSHION COVER,ASSY,SOL,4 WAY STRETCH,18"W,22"D,W/PVC	1
A1x	SOLSPP-NBLK1918	CUSHION COVER,ASSY,SOL,4 WAY STRETCH,19"W,18"D,W/PVC	1
A1y	SOLSPP-NBLK1919	CUSHION COVER,ASSY,SOL,4 WAY STRETCH,19"W,19"D,W/PVC	1
A1z	SOLSPP-NBLK1920	CUSHION COVER,ASSY,SOL,4 WAY STRETCH,19"W,20"D,W/PVC	1
A1a a	SOLSPP-NBLK2016	CUSHION COVER,ASSY,SOL,4 WAY STRETCH,20"W,16"D,W/PVC	1
A1b b	SOLSPP-NBLK2018	CUSHION COVER,ASSY,SOL,4 WAY STRETCH,20"W,18"D,W/PVC	1
A1c c	SOLSPP-NBLK2019	CUSHION COVER,ASSY,SOL,4 WAY STRETCH,20"W,19"D,W/PVC	1
A1d d	SOLSPP-NBLK2020	CUSHION COVER,ASSY,SOL,4 WAY STRETCH,20"W,20"D,W/PVC	1
A1e e	SOLSPP-NBLK2022	CUSHION COVER,ASSY,SOL,4 WAY STRETCH,20"W,22"D,W/PVC	1
A1ff	SOLSPP-NBLK2218	CUSHION COVER,ASSY,SOL,4 WAY STRETCH,22"W,18"D,W/PVC	1

**Stealth Solution  
4-way Stretch Outer Cushion Cover**

Size	Ref#	Size	Ref#	Size	Ref#
10"W x 10"D	A1a	16"W x 20"D	A1m	19"W x 19"D	A1y
10"W x 12"D	A1b	17"W x 16"D	A1n	19"W x 20"D	A1z
12"W x 10"D	A1c	17"W x 17"D	A1o	20"W x 16"D	A1aa
12"W x 12"D	A1d	17"W x 18"D	A1p	20"W x 18"D	A1bb
14"W x 12"D	A1e	17"W x 19"D	A1q	20"W x 19"D	A1cc
14"W x 14"D	A1f	17"W x 20"D	A1r	20"W x 20"D	A1dd
14"W x 16"D	A1g	18"W x 16"D	A1s	20"W x 22"D	A1ee
14"W x 18"D	A1h	18"W x 18"D	A1t	22"W x 18"D	A1ff
16"W x 16"D	A1i	18"W x 19"D	A1u	22"W x 20"D	A1gg
16"W x 17"D	A1j	18"W x 20"D	A1v	22"W x 22"D	A1hh
16"W x 18"D	A1k	18"W x 22"D	A1w	24"W x 20"D	A1ii
16"W x 19"D	A1l	19"W x 18"D	A1x	24"W x 22"D	A1jj






# Stealth Solution 4-way Stretch Outer Cushion Cover



Item	Part No.	Name	Qty
A1g	SOLSPP-NBLK2220	CUSHION COVER,ASSY,SOL,4 WAY STRETCH,22"W,20"D,W/PVC	1
A1h	SOLSPP-NBLK2222	CUSHION COVER,ASSY,SOL,4 WAY STRETCH,22"W,22"D,W/PVC	1
A1ii	SOLSPP-NBLK2420	CUSHION COVER,ASSY,SOL,4 WAY STRETCH,24"W,20"D,W/PVC	1
A1jj	SOLSPP-NBLK2422	CUSHION COVER,ASSY,SOL,4 WAY STRETCH,24"W,22"D,W/PVC	1

2. Outer Covers


Stealth Spectrum Foam 4-way Stretch Outer Cushion Cover							
Size	Ref#	Size	Ref#	Size	Ref#	Size	Ref#
16"W x 16"D	A1a	20"W x 16"D	A1m	24"W x 22"D	A1y	30"W x 20"D	A1kk
16"W x 18"D	A1b	20"W x 18"D	A1n	24"W x 24"D	A1z	30"W x 22"D	A1ll
16"W x 20"D	A1c	20"W x 20"D	A1o	26"W x 18"D	A1aa	30"W x 24"D	A1mm
17"W x 17"D	A1d	20"W x 22"D	A1p	26"W x 20"D	A1bb	30"W x 28"D	A1nn
17"W x 18"D	A1e	20"W x 24"D	A1q	26"W x 22"D	A1cc	32"W x 26"D	A1oo
17"W x 19"D	A1f	22"W x 18"D	A1r	26"W x 24"D	A1dd	32"W x 28"D	A1pp
18"W x 16"D	A1g	22"W x 20"D	A1s	28"W x 18"D	A1ee		
18"W x 18"D	A1h	22"W x 22"D	A1t	28"W x 20"D	A1ff		
18"W x 20"D	A1i	22"W x 24"D	A1u	28"W x 22"D	A1gg		
18"W x 22"D	A1j	24"W x 16"D	A1v	28"W x 24"D	A1hh		
19"W x 17"D	A1k	24"W x 18"D	A1w	28"W x 28"D	A1ii		
19"W x 19"D	A1l	24"W x 20"D	A1x	30"W x 18"D	A1jj		



Item	Part No.	Name	Qty
A1a	SPCFSPP-NBLK1616	CUSHION COVER,ASSY,SPCF,4 WAY STRETCH,16"W,16"D,W/PVC	1
A1b	SPCFSPP-NBLK1618	CUSHION COVER,ASSY,SPCF,4 WAY STRETCH,16"W,18"D,W/PVC	1
A1c	SPCFSPP-NBLK1620	CUSHION COVER,ASSY,SPCF,4 WAY STRETCH,16"W,20"D,W/PVC	1
A1d	SPCFSPP-NBLK1717	CUSHION COVER,ASSY,SPCF,4 WAY STRETCH,17"W,17"D,W/PVC	1
A1e	SPCFSPP-NBLK1718	CUSHION COVER,ASSY,SPCF,4 WAY STRETCH,17"W,18"D,W/PVC	1
A1f	SPCFSPP-NBLK1719	CUSHION COVER,ASSY,SPCF,4 WAY STRETCH,17"W,19"D,W/PVC	1
A1g	SPCFSPP-NBLK1816	CUSHION COVER,ASSY,SPCF,4 WAY STRETCH,18"W,16"D,W/PVC	1
A1h	SPCFSPP-NBLK1818	CUSHION COVER,ASSY,SPCF,4 WAY STRETCH,18"W,18"D,W/PVC	1
A1i	SPCFSPP-NBLK1820	CUSHION COVER,ASSY,SPCF,4 WAY STRETCH,18"W,20"D,W/PVC	1
A1j	SPCFSPP-NBLK1822	CUSHION COVER,ASSY,SPCF,4 WAY STRETCH,18"W,22"D,W/PVC	1
A1k	SPCFSPP-NBLK1917	CUSHION COVER,ASSY,SPCF,4 WAY STRETCH,19"W,17"D,W/PVC	1
A1l	SPCFSPP-NBLK1919	CUSHION COVER,ASSY,SPCF,4 WAY STRETCH,19"W,19"D,W/PVC	1
A1m	SPCFSPP-NBLK2016	CUSHION COVER,ASSY,SPCF,4 WAY STRETCH,20"W,16"D,W/PVC	1
A1n	SPCFSPP-NBLK2018	CUSHION COVER,ASSY,SPCF,4 WAY STRETCH,20"W,18"D,W/PVC	1
A1o	SPCFSPP-NBLK2020	CUSHION COVER,ASSY,SPCF,4 WAY STRETCH,20"W,20"D,W/PVC	1
A1p	SPCFSPP-NBLK2022	CUSHION COVER,ASSY,SPCF,4 WAY STRETCH,20"W,22"D,W/PVC	1
A1q	SPCFSPP-NBLK2024	CUSHION COVER,ASSY,SPCF,4 WAY STRETCH,20"W,24"D,W/PVC	1
A1r	SPCFSPP-NBLK2218	CUSHION COVER,ASSY,SPCF,4 WAY STRETCH,22"W,18"D,W/PVC	1
A1s	SPCFSPP-NBLK2220	CUSHION COVER,ASSY,SPCF,4 WAY STRETCH,22"W,20"D,W/PVC	1
A1t	SPCFSPP-NBLK2222	CUSHION COVER,ASSY,SPCF,4 WAY STRETCH,22"W,22"D,W/PVC	1
A1u	SPCFSPP-NBLK2224	CUSHION COVER,ASSY,SPCF,4 WAY STRETCH,22"W,24"D,W/PVC	1
A1v	SPCFSPP-NBLK2416	CUSHION COVER,ASSY,SPCF,4 WAY STRETCH,24"W,16"D,W/PVC	1
A1w	SPCFSPP-NBLK2418	CUSHION COVER,ASSY,SPCF,4 WAY STRETCH,24"W,18"D,W/PVC	1
A1x	SPCFSPP-NBLK2420	CUSHION COVER,ASSY,SPCF,4 WAY STRETCH,24"W,20"D,W/PVC	1

### Stealth Spectrum Foam 4-way Stretch Outer Cushion Cover

Size	Ref#	Size	Ref#	Size	Ref#	Size	Ref#
16"W x 16"D	A1a	20"W x 16"D	A1m	24"W x 22"D	A1y	30"W x 20"D	A1kk
16"W x 18"D	A1b	20"W x 18"D	A1n	24"W x 24"D	A1z	30"W x 22"D	A1ll
16"W x 20"D	A1c	20"W x 20"D	A1o	26"W x 18"D	A1aa	30"W x 24"D	A1mm
17"W x 17"D	A1d	20"W x 22"D	A1p	26"W x 20"D	A1bb	30"W x 28"D	A1nn
17"W x 18"D	A1e	20"W x 24"D	A1q	26"W x 22"D	A1cc	32"W x 26"D	A1oo
17"W x 19"D	A1f	22"W x 18"D	A1r	26"W x 24"D	A1dd	32"W x 28"D	A1pp
18"W x 16"D	A1g	22"W x 20"D	A1s	28"W x 18"D	A1ee		
18"W x 18"D	A1h	22"W x 22"D	A1t	28"W x 20"D	A1ff		
18"W x 20"D	A1i	22"W x 24"D	A1u	28"W x 22"D	A1gg		
18"W x 22"D	A1j	24"W x 16"D	A1v	28"W x 24"D	A1hh		
19"W x 17"D	A1k	24"W x 18"D	A1w	28"W x 28"D	A1ii		
19"W x 19"D	A1l	24"W x 20"D	A1x	30"W x 18"D	A1jj		




Item	Part No.	Name	Qty
A1y	SPCFSP- NBLK2422	CUSHION COVER,ASSY,SPCF,4 WAY STRETCH,24"W,22"D,W/PVC	1
A1z	SPCFSP- NBLK2424	CUSHION COVER,ASSY,SPCF,4 WAY STRETCH,24"W,24"D,W/PVC	1
A1a	SPCFSP- NBLK2618	CUSHION COVER,ASSY,SPCF,4 WAY STRETCH,26"W,18"D,W/PVC	1
A1b	SPCFSP- NBLK2620	CUSHION COVER,ASSY,SPCF,4 WAY STRETCH,26"W,20"D,W/PVC	1
A1c	SPCFSP- NBLK2622	CUSHION COVER,ASSY,SPCF,4 WAY STRETCH,26"W,22"D,W/PVC	1
A1d	SPCFSP- NBLK2624	CUSHION COVER,ASSY,SPCF,4 WAY STRETCH,26"W,24"D,W/PVC	1
A1e	SPCFSP- NBLK2818	CUSHION COVER,ASSY,SPCF,4 WAY STRETCH,28"W,18"D,W/PVC	1
A1ff	SPCFSP- NBLK2820	CUSHION COVER,ASSY,SPCF,4 WAY STRETCH,28"W,20"D,W/PVC	1
A1g	SPCFSP- NBLK2822	CUSHION COVER,ASSY,SPCF,4 WAY STRETCH,28"W,22"D,W/PVC	1
A1h	SPCFSP- NBLK2824	CUSHION COVER,ASSY,SPCF,4 WAY STRETCH,28"W,24"D,W/PVC	1
A1ii	SPCFSP- NBLK2828	CUSHION COVER,ASSY,SPCF,4 WAY STRETCH,28"W,28"D,W/PVC	1
A1jj	SPCFSP- NBLK3018	CUSHION COVER,ASSY,SPCF,4 WAY STRETCH,30"W,18"D,W/PVC	1
A1k	SPCFSP- NBLK3020	CUSHION COVER,ASSY,SPCF,4 WAY STRETCH,30"W,20"D,W/PVC	1
A1ll	SPCFSP- NBLK3022	CUSHION COVER,ASSY,SPCF,4 WAY STRETCH,30"W,22"D,W/PVC	1
A1m	SPCFSP- NBLK3024	CUSHION COVER,ASSY,SPCF,4 WAY STRETCH,30"W,24"D,W/PVC	1
A1n	SPCFSP- NBLK3028	CUSHION COVER,ASSY,SPCF,4 WAY STRETCH,30"W,28"D,W/PVC	1
A1o	SPCFSP- NBLK3226	CUSHION COVER,ASSY,SPCF,4 WAY STRETCH,32"W,26"D,W/PVC	1
A1p	SPCFSP- NBLK3228	CUSHION COVER,ASSY,SPCF,4 WAY STRETCH,32"W,28"D,W/PVC	1

2. Outer Covers

**Stealth Spectrum Gel  
Outer Cushion Cover**


Size	Ref#	Size	Ref#	Size	Ref#
14"W x 18"D	A1a	18"W x 20"D	A1n	22"W x 16"D	A1aa
16"W x 16"D	A1b	18"W x 22"D	A1o	22"W x 18"D	A1bb
16"W x 18"D	A1c	18"W x 24"D	A1p	22"W x 19"D	A1cc
16"W x 19"D	A1d	19"W x 17"D	A1q	22"W x 20"D	A1dd
16"W x 20"D	A1e	19"W x 18"D	A1r	22"W x 22"D	A1ee
17"W x 16"D	A1f	19"W x 19"D	A1s	22"W x 24"D	A1ff
17"W x 17"D	A1g	19"W x 20"D	A1t	24"W x 16"D	A1gg
17"W x 18"D	A1h	20"W x 16"D	A1u	24"W x 18"D	A1hh
17"W x 19"D	A1i	20"W x 18"D	A1v	24"W x 20"D	A1ii
17"W x 20"D	A1j	20"W x 19"D	A1w	24"W x 22"D	A1jj
18"W x 16"D	A1k	20"W x 20"D	A1x	24"W x 24"D	A1kk
18"W x 17"D	A1l	20"W x 22"D	A1y		
18"W x 18"D	A1m	20"W x 24"D	A1z		



Item	Part No.	Name	Qty
A1a	SPCGSPP-NBLK1418	CUSHION COVER,ASSY,SPCG,4 WAY STRETCH,14"W,18"D,W/PVC	1
A1b	SPCGSPP-NBLK1616	CUSHION COVER,ASSY,SPCG,4 WAY STRETCH,16"W,16"D,W/PVC	1
A1c	SPCGSPP-NBLK1618	CUSHION COVER,ASSY,SPCG,4 WAY STRETCH,16"W,18"D,W/PVC	1
A1d	SPCGSPP-NBLK1619	CUSHION COVER,ASSY,SPCG,4 WAY STRETCH,16"W,19"D,W/PVC	1
A1e	SPCGSPP-NBLK1620	CUSHION COVER,ASSY,SPCG,4 WAY STRETCH,16"W,20"D,W/PVC	1
A1f	SPCGSPP-NBLK1716	CUSHION COVER,ASSY,SPCG,4 WAY STRETCH,17"W,16"D,W/PVC	1
A1g	SPCGSPP-NBLK1717	CUSHION COVER,ASSY,SPCG,4 WAY STRETCH,17"W,17"D,W/PVC	1
A1h	SPCGSPP-NBLK1718	CUSHION COVER,ASSY,SPCG,4 WAY STRETCH,17"W,18"D,W/PVC	1
A1i	SPCGSPP-NBLK1719	CUSHION COVER,ASSY,SPCG,4 WAY STRETCH,17"W,19"D,W/PVC	1
A1j	SPCGSPP-NBLK1720	CUSHION COVER,ASSY,SPCG,4 WAY STRETCH,17"W,20"D,W/PVC	1
A1k	SPCGSPP-NBLK1816	CUSHION COVER,ASSY,SPCG,4 WAY STRETCH,18"W,16"D,W/PVC	1
A1l	SPCGSPP-NBLK1817	CUSHION COVER,ASSY,SPCG,4 WAY STRETCH,18"W,17"D,W/PVC	1
A1m	SPCGSPP-NBLK1818	CUSHION COVER,ASSY,SPCG,4 WAY STRETCH,18"W,18"D,W/PVC	1
A1n	SPCGSPP-NBLK1820	CUSHION COVER,ASSY,SPCG,4 WAY STRETCH,18"W,20"D,W/PVC	1
A1o	SPCGSPP-NBLK1822	CUSHION COVER,ASSY,SPCG,4 WAY STRETCH,18"W,22"D,W/PVC	1
A1p	SPCGSPP-NBLK1824	CUSHION COVER,ASSY,SPCG,4 WAY STRETCH,18"W,24"D,W/PVC	1
A1q	SPCGSPP-NBLK1917	CUSHION COVER,ASSY,SPCG,4 WAY STRETCH,19"W,17"D,W/PVC	1
A1r	SPCGSPP-NBLK1918	CUSHION COVER,ASSY,SPCG,4 WAY STRETCH,19"W,18"D,W/PVC	1
A1s	SPCGSPP-NBLK1919	CUSHION COVER,ASSY,SPCG,4 WAY STRETCH,19"W,19"D,W/PVC	1
A1t	SPCGSPP-NBLK1920	CUSHION COVER,ASSY,SPCG,4 WAY STRETCH,19"W,20"D,W/PVC	1
A1u	SPCGSPP-NBLK2016	CUSHION COVER,ASSY,SPCG,4 WAY STRETCH,20"W,16"D,W/PVC	1
A1v	SPCGSPP-NBLK2018	CUSHION COVER,ASSY,SPCG,4 WAY STRETCH,20"W,18"D,W/PVC	1
A1w	SPCGSPP-NBLK2019	CUSHION COVER,ASSY,SPCG,4 WAY STRETCH,20"W,19"D,W/PVC	1
A1x	SPCGSPP-NBLK2020	CUSHION COVER,ASSY,SPCG,4 WAY STRETCH,20"W,20"D,W/PVC	1

### Stealth Spectrum Gel Outer Cushion Cover

Size	Ref#	Size	Ref#	Size	Ref#
14"W x 18"D	A1a	18"W x 20"D	A1n	22"W x 16"D	A1aa
16"W x 16"D	A1b	18"W x 22"D	A1o	22"W x 18"D	A1bb
16"W x 18"D	A1c	18"W x 24"D	A1p	22"W x 19"D	A1cc
16"W x 19"D	A1d	19"W x 17"D	A1q	22"W x 20"D	A1dd
16"W x 20"D	A1e	19"W x 18"D	A1r	22"W x 22"D	A1ee
17"W x 16"D	A1f	19"W x 19"D	A1s	22"W x 24"D	A1ff
17"W x 17"D	A1g	19"W x 20"D	A1t	24"W x 16"D	A1gg
17"W x 18"D	A1h	20"W x 16"D	A1u	24"W x 18"D	A1hh
17"W x 19"D	A1i	20"W x 18"D	A1v	24"W x 20"D	A1ii
17"W x 20"D	A1j	20"W x 19"D	A1w	24"W x 22"D	A1jj
18"W x 16"D	A1k	20"W x 20"D	A1x	24"W x 24"D	A1kk
18"W x 17"D	A1l	20"W x 22"D	A1y		
18"W x 18"D	A1m	20"W x 24"D	A1z		






Item	Part No.	Name	Qty
A1y	SPCGSPP-NBLK2022	CUSHION COVER,ASSY,SPCG,4 WAY STRETCH,20"W,22"D,W/PVC	1
A1z	SPCGSPP-NBLK2024	CUSHION COVER,ASSY,SPCG,4 WAY STRETCH,20"W,24"D,W/PVC	1
A1a	SPCGSPP-NBLK2216	CUSHION COVER,ASSY,SPCG,4 WAY STRETCH,22"W,16"D,W/PVC	1
A1b	SPCGSPP-NBLK2218	CUSHION COVER,ASSY,SPCG,4 WAY STRETCH,22"W,18"D,W/PVC	1
A1c	SPCGSPP-NBLK2219	CUSHION COVER,ASSY,SPCG,4 WAY STRETCH,22"W,19"D,W/PVC	1
A1d	SPCGSPP-NBLK2220	CUSHION COVER,ASSY,SPCG,4 WAY STRETCH,22"W,20"D,W/PVC	1
A1e	SPCGSPP-NBLK2222	CUSHION COVER,ASSY,SPCG,4 WAY STRETCH,22"W,22"D,W/PVC	1
A1ff	SPCGSPP-NBLK2224	CUSHION COVER,ASSY,SPCG,4 WAY STRETCH,22"W,24"D,W/PVC	1
A1g	SPCGSPP-NBLK2416	CUSHION COVER,ASSY,SPCG,4 WAY STRETCH,24"W,16"D,W/PVC	1
A1h	SPCGSPP-NBLK2418	CUSHION COVER,ASSY,SPCG,4 WAY STRETCH,24"W,18"D,W/PVC	1
A1ii	SPCGSPP-NBLK2420	CUSHION COVER,ASSY,SPCG,4 WAY STRETCH,24"W,20"D,W/PVC	1
A1jj	SPCGSPP-NBLK2422	CUSHION COVER,ASSY,SPCG,4 WAY STRETCH,24"W,22"D,W/PVC	1
A1k	SPCGSPP-NBLK2424	CUSHION COVER,ASSY,SPCG,4 WAY STRETCH,24"W,24"D,W/PVC	1

2. Outer Covers

Stealth TRU-Comfort 2 Outer Cushion Cover			
Size	Ref#	Size	Ref#
16"W x 14"D	A1a	20"W x 18"D	A1j
16"W x 16"D	A1b	20"W x 20"D	A1k
16"W x 18"D	A1c	20"W x 22"D	A1l
16"W x 20"D	A1d	20"W x 24"D	A1m
18"W x 16"D	A1e	22"W x 18"D	A1n
18"W x 18"D	A1f	22"W x 20"D	A1o
18"W x 20"D	A1g	22"W x 22"D	A1p
18"W x 22"D	A1h	22"W x 24"D	A1q
20"W x 16"D	A1i		




Item	Part No.	Name	Qty
A1a	TC2SPP-NBLK1614	CUSHION COVER,ASSY,TC2,4 WAY STRETCH,16"W,14"D,W/PVC	1
A1b	TC2SPP-NBLK1616	CUSHION COVER,ASSY,TC2,4 WAY STRETCH,16"W,16"D,W/PVC	1
A1c	TC2SPP-NBLK1618	CUSHION COVER,ASSY,TC2,4 WAY STRETCH,16"W,18"D,W/PVC	1
A1d	TC2SPP-NBLK1620	CUSHION COVER,ASSY,TC2,4 WAY STRETCH,16"W,20"D,W/PVC	1
A1e	TC2SPP-NBLK1816	CUSHION COVER,ASSY,TC2,4 WAY STRETCH,18"W,16"D,W/PVC	1
A1f	TC2SPP-NBLK1818	CUSHION COVER,ASSY,TC2,4 WAY STRETCH,18"W,18"D,W/PVC	1
A1g	TC2SPP-NBLK1820	CUSHION COVER,ASSY,TC2,4 WAY STRETCH,18"W,20"D,W/PVC	1
A1h	TC2SPP-NBLK1822	CUSHION COVER,ASSY,TC2,4 WAY STRETCH,18"W,22"D,W/PVC	1
A1i	TC2SPP-NBLK2016	CUSHION COVER,ASSY,TC2,4 WAY STRETCH,20"W,16"D,W/PVC	1
A1j	TC2SPP-NBLK2018	CUSHION COVER,ASSY,TC2,4 WAY STRETCH,20"W,18"D,W/PVC	1
A1k	TC2SPP-NBLK2020	CUSHION COVER,ASSY,TC2,4 WAY STRETCH,20"W,20"D,W/PVC	1
A1l	TC2SPP-NBLK2022	CUSHION COVER,ASSY,TC2,4 WAY STRETCH,20"W,22"D,W/PVC	1
A1m	TC2SPP-NBLK2024	CUSHION COVER,ASSY,TC2,4 WAY STRETCH,20"W,24"D,W/PVC	1
A1n	TC2SPP-NBLK2218	CUSHION COVER,ASSY,TC2,4 WAY STRETCH,22"W,18"D,W/PVC	1
A1o	TC2SPP-NBLK2220	CUSHION COVER,ASSY,TC2,4 WAY STRETCH,22"W,20"D,W/PVC	1
A1p	TC2SPP-NBLK2222	CUSHION COVER,ASSY,TC2,4 WAY STRETCH,22"W,22"D,W/PVC	1
A1q	TC2SPP-NBLK2224	CUSHION COVER,ASSY,TC2,4 WAY STRETCH,22"W,24"D,W/PVC	1

### 3. Inner Covers

**Urethane Inner Liner - 10"W - 19"W  
Used with Solution, Spectrum Gel,  
and Spectrum Foam Cushions**

Size	Ref#	Size	Ref#	Size	Ref#
10"W x 10"D	A1a	16"W x 18"D	A1k	18"W x 18"D	A1u
10"W x 12"D	A1b	16"W x 19"D	A1l	18"W x 19"D	A1v
12"W x 10"D	A1c	16"W x 20"D	A1m	18"W x 20"D	A1w
12"W x 12"D	A1d	17"W x 16"D	A1n	18"W x 22"D	A1x
14"W x 12"D	A1e	17"W x 17"D	A1o	18"W x 24"D	A1y
14"W x 14"D	A1f	17"W x 18"D	A1p	19"W x 17"D	A1z
14"W x 16"D	A1g	17"W x 19"D	A1q	19"W x 18"D	A1aa
14"W x 18"D	A1h	17"W x 20"D	A1r	19"W x 19"D	A1bb
16"W x 16"D	A1i	18"W x 16"D	A1s	19"W x 20"D	A1cc
16"W x 17"D	A1j	18"W x 17"D	A1t		

\*The urethane liner is a standard component for the cushions listed and is included with the top level cushion assemblies.




Item	Part No.	Name	Qty
A1a	URL1709070-1010	LINER,URETHANE,SOL,SOL1,SPCF,SPCG,10"W,10"D	1
A1b	URL1709070-1012	LINER,URETHANE,SOL,SOL1,SPCF,SPCG,10"W,12"D	1
A1c	URL1709070-1210	LINER,URETHANE,SOL,SOL1,SPCF,SPCG,12"W,10"D	1
A1d	URL1709070-1212	LINER,URETHANE,SOL,SOL1,SPCF,SPCG,12"W,12"D	1
A1e	URL1709070-1412	LINER,URETHANE,SOL,SOL1,SPCF,SPCG,14"W,12"D	1
A1f	URL1709070-1414	LINER,URETHANE,SOL,SOL1,SPCF,SPCG,14"W,14"D	1
A1g	URL1709070-1416	LINER,URETHANE,SOL,SOL1,SPCF,SPCG,14"W,16"D	1
A1h	URL1709070-1418	LINER,URETHANE,SOL,SOL1,SPCF,SPCG,14"W,18"D	1
A1i	URL1709070-1616	LINER,URETHANE,SOL,SOL1,SPCF,SPCG,16"W,16"D	1
A1j	URL1709070-1617	LINER,URETHANE,SOL,SOL1,SPCF,SPCG,16"W,17"D	1
A1k	URL1709070-1618	LINER,URETHANE,SOL,SOL1,SPCF,SPCG,16"W,18"D	1
A1l	URL1709070-1619	LINER,URETHANE,SOL,SOL1,SPCF,SPCG,16"W,19"D	1
A1m	URL1709070-1620	LINER,URETHANE,SOL,SOL1,SPCF,SPCG,16"W,20"D	1
A1n	URL1709070-1716	LINER,URETHANE,SOL,SOL1,SPCF,SPCG,17"W,16"D	1
A1o	URL1709070-1717	LINER,URETHANE,SOL,SOL1,SPCF,SPCG,17"W,17"D	1
A1p	URL1709070-1718	LINER,URETHANE,SOL,SOL1,SPCF,SPCG,17"W,18"D	1
A1q	URL1709070-1719	LINER,URETHANE,SOL,SOL1,SPCF,SPCG,17"W,19"D	1
A1r	URL1709070-1720	LINER,URETHANE,SOL,SOL1,SPCF,SPCG,17"W,20"D	1
A1s	URL1709070-1816	LINER,URETHANE,SOL,SOL1,SPCF,SPCG,18"W,16"D	1
A1t	URL1709070-1817	LINER,URETHANE,SOL,SOL1,SPCF,SPCG,18"W,17"D	1
A1u	URL1709070-1818	LINER,URETHANE,SOL,SOL1,SPCF,SPCG,18"W,18"D	1
A1v	URL1709070-1819	LINER,URETHANE,SOL,SOL1,SPCF,SPCG,18"W,19"D	1
A1w	URL1709070-1820	LINER,URETHANE,SOL,SOL1,SPCF,SPCG,18"W,20"D	1
A1x	URL1709070-1822	LINER,URETHANE,SOL,SOL1,SPCF,SPCG,18"W,22"D	1
A1y	URL1709070-1824	LINER,URETHANE,SOL,SOL1,SPCF,SPCG,18"W,24"D	1
A1z	URL1709070-1917	LINER,URETHANE,SOL,SOL1,SPCF,SPCG,19"W,17"D	1
A1a a	URL1709070-1918	LINER,URETHANE,SOL,SOL1,SPCF,SPCG,19"W,18"D	1
A1b b	URL1709070-1919	LINER,URETHANE,SOL,SOL1,SPCF,SPCG,19"W,19"D	1
A1c c	URL1709070-1920	LINER,URETHANE,SOL,SOL1,SPCF,SPCG,19"W,20"D	1

### 3. Inner Covers

**Urethane Inner Liner - 20"W - 32"W  
Used with Solution, Spectrum Gel,  
and Spectrum Foam Cushions**

Size	Ref#	Size	Ref#	Size	Ref#
20"W x 16"D	A1a	22"W x 24"D	A1l	28"W x 20"D	A1w
20"W x 18"D	A1b	24"W x 16"D	A1m	28"W x 22"D	A1x
20"W x 19"D	A1c	24"W x 18"D	A1n	28"W x 24"D	A1y
20"W x 20"D	A1d	24"W x 20"D	A1o	28"W x 28"D	A1z
20"W x 22"D	A1e	24"W x 22"D	A1p	30"W x 18"D	A1aa
20"W x 24"D	A1f	24"W x 24"D	A1q	30"W x 20"D	A1bb
22"W x 16"D	A1g	26"W x 18"D	A1r	30"W x 22"D	A1cc
22"W x 18"D	A1h	26"W x 20"D	A1s	30"W x 24"D	A1dd
22"W x 19"D	A1i	26"W x 22"D	A1t	30"W x 28"D	A1ee
22"W x 20"D	A1j	26"W x 24"D	A1u	32"W x 26"D	A1ff
22"W x 22"D	A1k	28"W x 18"D	A1v	32"W x 28"D	A1gg

\*The urethane liner is a standard component for the cushions listed and is included with the top level cushion assemblies.




Ite m	Part No.	Name	Qty
A1a	URL1709070-2016	LINER,URETHANE,SOL,SOL1,SPCF,SPCG,20"W,16"D	1
A1b	URL1709070-2018	LINER,URETHANE,SOL,SOL1,SPCF,SPCG,20"W,18"D	1
A1c	URL1709070-2019	LINER,URETHANE,SOL,SOL1,SPCF,SPCG,20"W,19"D	1
A1d	URL1709070-2020	LINER,URETHANE,SOL,SOL1,SPCF,SPCG,20"W,20"D	1
A1e	URL1709070-2022	LINER,URETHANE,SOL,SOL1,SPCF,SPCG,20"W,22"D	1
A1f	URL1709070-2024	LINER,URETHANE,SOL,SOL1,SPCF,SPCG,20"W,24"D	1
A1g	URL1709070-2216	LINER,URETHANE,SOL,SOL1,SPCF,SPCG,22"W,16"D	1
A1h	URL1709070-2218	LINER,URETHANE,SOL,SOL1,SPCF,SPCG,22"W,18"D	1
A1i	URL1709070-2219	LINER,URETHANE,SOL,SOL1,SPCF,SPCG,22"W,19"D	1
A1j	URL1709070-2220	LINER,URETHANE,SOL,SOL1,SPCF,SPCG,22"W,20"D	1
A1k	URL1709070-2222	LINER,URETHANE,SOL,SOL1,SPCF,SPCG,22"W,22"D	1
A1l	URL1709070-2224	LINER,URETHANE,SOL,SOL1,SPCF,SPCG,22"W,24"D	1
A1m	URL1709070-2416	LINER,URETHANE,SOL,SOL1,SPCF,SPCG,24"W,16"D	1
A1n	URL1709070-2418	LINER,URETHANE,SOL,SOL1,SPCF,SPCG,24"W,18"D	1
A1o	URL1709070-2420	LINER,URETHANE,SOL,SOL1,SPCF,SPCG,24"W,20"D	1
A1p	URL1709070-2422	LINER,URETHANE,SOL,SOL1,SPCF,SPCG,24"W,22"D	1
A1q	URL1709070-2424	LINER,URETHANE,SOL,SOL1,SPCF,SPCG,24"W,24"D	1
A1r	URL1709070-2618	LINER,URETHANE,SOL,SOL1,SPCF,SPCG,26"W,18"D	1
A1s	URL1709070-2620	LINER,URETHANE,SOL,SOL1,SPCF,SPCG,26"W,20"D	1
A1t	URL1709070-2622	LINER,URETHANE,SOL,SOL1,SPCF,SPCG,26"W,22"D	1
A1u	URL1709070-2624	LINER,URETHANE,SOL,SOL1,SPCF,SPCG,26"W,24"D	1
A1v	URL1709070-2818	LINER,URETHANE,SOL,SOL1,SPCF,SPCG,28"W,18"D	1
A1w	URL1709070-2820	LINER,URETHANE,SOL,SOL1,SPCF,SPCG,28"W,20"D	1
A1x	URL1709070-2822	LINER,URETHANE,SOL,SOL1,SPCF,SPCG,28"W,22"D	1
A1y	URL1709070-2824	LINER,URETHANE,SOL,SOL1,SPCF,SPCG,28"W,24"D	1
A1z	URL1709070-2828	LINER,URETHANE,SOL,SOL1,SPCF,SPCG,28"W,28"D	1
A1a a	URL1709070-3018	LINER,URETHANE,SOL,SOL1,SPCF,SPCG,30"W,18"D	1
A1b b	URL1709070-3020	LINER,URETHANE,SOL,SOL1,SPCF,SPCG,30"W,20"D	1
A1c c	URL1709070-3022	LINER,URETHANE,SOL,SOL1,SPCF,SPCG,30"W,22"D	1
A1d d	URL1709070-3024	LINER,URETHANE,SOL,SOL1,SPCF,SPCG,30"W,24"D	1
A1e e	URL1709070-3028	LINER,URETHANE,SOL,SOL1,SPCF,SPCG,30"W,28"D	1
A1ff	URL1709070-3226	LINER,URETHANE,SOL,SOL1,SPCF,SPCG,32"W,26"D	1

**Urethane Inner Liner - 20"W - 32"W  
Used with Solution, Spectrum Gel,  
and Spectrum Foam Cushions**

Size	Ref#	Size	Ref#	Size	Ref#
20"W x 16"D	A1a	22"W x 24"D	A1l	28"W x 20"D	A1w
20"W x 18"D	A1b	24"W x 16"D	A1m	28"W x 22"D	A1x
20"W x 19"D	A1c	24"W x 18"D	A1n	28"W x 24"D	A1y
20"W x 20"D	A1d	24"W x 20"D	A1o	28"W x 28"D	A1z
20"W x 22"D	A1e	24"W x 22"D	A1p	30"W x 18"D	A1aa
20"W x 24"D	A1f	24"W x 24"D	A1q	30"W x 20"D	A1bb
22"W x 16"D	A1g	26"W x 18"D	A1r	30"W x 22"D	A1cc
22"W x 18"D	A1h	26"W x 20"D	A1s	30"W x 24"D	A1dd
22"W x 19"D	A1i	26"W x 22"D	A1t	30"W x 28"D	A1ee
22"W x 20"D	A1j	26"W x 24"D	A1u	32"W x 26"D	A1ff
22"W x 22"D	A1k	28"W x 18"D	A1v	32"W x 28"D	A1gg



\*The urethane liner is a standard component for the cushions listed and is included with the top level cushion assemblies.





Item	Part No.	Name	Qty
A1gg	URL1709070-3228	LINER,URETHANE,SOL,SOL1,SPCF,SPCG,32"W,28"D	1

### 3. Inner Covers



Size	Ref#	Size	Ref#
16"W x 16"D	A1	18"W x 20"D	A5
16"W x 18"D	A2	20"W x 18"D	A6
18"W x 16"D	A3	20"W x 20"D	A7
18"W x 18"D	A4		

\*Please note: The polyurethane inner liner is an optional component and not included in the top level cushion assemblies.


Item	Part No.	Name	Qty
A1	PUL1709080-1616	LINER,POLYURETHANE,TC2,ESS,16"W,16"D	1
A2	PUL1709080-1618	LINER,POLYURETHANE,TC2,ESS,16"W,18"D	1
A3	PUL1709080-1816	LINER,POLYURETHANE,TC2,ESS,18"W,16"D	1
A4	PUL1709080-1818	LINER,POLYURETHANE,TC2,ESS,18"W,18"D	1
A5	PUL1709080-1820	LINER,POLYURETHANE,TC2,ESS,18"W,20"D	1
A6	PUL1709080-2018	LINER,POLYURETHANE,TC2,ESS,20"W,18"D	1
A7	PUL1709080-2020	LINER,POLYURETHANE,TC2,ESS,20"W,20"D	1

### 3. Inner Covers

**Polyurethane Inner Liner - 10"W - 19"W  
Used with Stealth Solution, Spectrum Gel, and  
Spectrum Foam Cushions**

Size	Ref#	Size	Ref#	Size	Ref#
10"W x 10"D	A1a	16"W x 18"D	A1k	18"W x 18"D	A1u
10"W x 12"D	A1b	16"W x 19"D	A1l	18"W x 19"D	A1v
12"W x 10"D	A1c	16"W x 20"D	A1m	18"W x 20"D	A1w
12"W x 12"D	A1d	17"W x 16"D	A1n	18"W x 22"D	A1x
14"W x 12"D	A1e	17"W x 17"D	A1o	18"W x 24"D	A1y
14"W x 14"D	A1f	17"W x 18"D	A1p	19"W x 17"D	A1z
14"W x 16"D	A1g	17"W x 19"D	A1q	19"W x 18"D	A1aa
14"W x 18"D	A1h	17"W x 20"D	A1r	19"W x 19"D	A1bb
16"W x 16"D	A1i	18"W x 16"D	A1s	19"W x 20"D	A1cc
16"W x 17"D	A1j	18"W x 17"D	A1t		

\*Please note: The polyurethane inner liner is an optional component and not included in the top level cushion assemblies.




Ite m	Part No.	Name	Qty
A1a	PUL1709030-1010	LINER,POLYURETHANE,SOL,SOL1,SPCF,SPCG,10"W,10"D	1
A1b	PUL1709030-1012	LINER,POLYURETHANE,SOL,SOL1,SPCF,SPCG,10"W,12"D	1
A1c	PUL1709030-1210	LINER,POLYURETHANE,SOL,SOL1,SPCF,SPCG,12"W,10"D	1
A1d	PUL1709030-1212	LINER,POLYURETHANE,SOL,SOL1,SPCF,SPCG,12"W,12"D	1
A1e	PUL1709030-1412	LINER,POLYURETHANE,SOL,SOL1,SPCF,SPCG,14"W,12"D	1
A1f	PUL1709030-1414	LINER,POLYURETHANE,SOL,SOL1,SPCF,SPCG,14"W,14"D	1
A1g	PUL1709030-1416	LINER,POLYURETHANE,SOL,SOL1,SPCF,SPCG,14"W,16"D	1
A1h	PUL1709030-1418	LINER,POLYURETHANE,SOL,SOL1,SPCF,SPCG,14"W,18"D	1
A1i	PUL1709030-1616	LINER,POLYURETHANE,SOL,SOL1,SPCF,SPCG,16"W,16"D	1
A1j	PUL1709030-1617	LINER,POLYURETHANE,SOL,SOL1,SPCF,SPCG,16"W,17"D	1
A1k	PUL1709030-1618	LINER,POLYURETHANE,SOL,SOL1,SPCF,SPCG,16"W,18"D	1
A1l	PUL1709030-1619	LINER,POLYURETHANE,SOL,SOL1,SPCF,SPCG,16"W,19"D	1
A1m	PUL1709030-1620	LINER,POLYURETHANE,SOL,SOL1,SPCF,SPCG,16"W,20"D	1
A1n	PUL1709030-1716	LINER,POLYURETHANE,SOL,SOL1,SPCF,SPCG,17"W,16"D	1
A1o	PUL1709030-1717	LINER,POLYURETHANE,SOL,SOL1,SPCF,SPCG,17"W,17"D	1
A1p	PUL1709030-1718	LINER,POLYURETHANE,SOL,SOL1,SPCF,SPCG,17"W,18"D	1
A1q	PUL1709030-1719	LINER,POLYURETHANE,SOL,SOL1,SPCF,SPCG,17"W,19"D	1
A1r	PUL1709030-1720	LINER,POLYURETHANE,SOL,SOL1,SPCF,SPCG,17"W,20"D	1
A1s	PUL1709030-1816	LINER,POLYURETHANE,SOL,SOL1,SPCF,SPCG,18"W,16"D	1
A1t	PUL1709030-1817	LINER,POLYURETHANE,SOL,SOL1,SPCF,SPCG,18"W,17"D	1
A1u	PUL1709030-1818	LINER,POLYURETHANE,SOL,SOL1,SPCF,SPCG,18"W,18"D	1
A1v	PUL1709030-1819	LINER,POLYURETHANE,SOL,SOL1,SPCF,SPCG,18"W,19"D	1
A1w	PUL1709030-1820	LINER,POLYURETHANE,SOL,SOL1,SPCF,SPCG,18"W,20"D	1
A1x	PUL1709030-1822	LINER,POLYURETHANE,SOL,SOL1,SPCF,SPCG,18"W,22"D	1
A1y	PUL1709030-1824	LINER,POLYURETHANE,SOL,SOL1,SPCF,SPCG,18"W,24"D	1
A1z	PUL1709030-1917	LINER,POLYURETHANE,SOL,SOL1,SPCF,SPCG,19"W,17"D	1
A1a a	PUL1709030-1918	LINER,POLYURETHANE,SOL,SOL1,SPCF,SPCG,19"W,18"D	1
A1b b	PUL1709030-1919	LINER,POLYURETHANE,SOL,SOL1,SPCF,SPCG,19"W,19"D	1
A1c c	PUL1709030-1920	LINER,POLYURETHANE,SOL,SOL1,SPCF,SPCG,19"W,20"D	1

### 3. Inner Covers

**Polyurethane Inner Liner - 20"W - 32"W  
Used with Solution, Spectrum Gel,  
and Spectrum Foam Cushions**

Size	Ref#	Size	Ref#	Size	Ref#
20"W x 16"D	A1a	22"W x 24"D	A1l	28"W x 20"D	A1w
20"W x 18"D	A1b	24"W x 16"D	A1m	28"W x 22"D	A1x
20"W x 19"D	A1c	24"W x 18"D	A1n	28"W x 24"D	A1y
20"W x 20"D	A1d	24"W x 20"D	A1o	28"W x 28"D	A1z
20"W x 22"D	A1e	24"W x 22"D	A1p	30"W x 18"D	A1aa
20"W x 24"D	A1f	24"W x 24"D	A1q	30"W x 20"D	A1bb
22"W x 16"D	A1g	26"W x 18"D	A1r	30"W x 22"D	A1cc
22"W x 18"D	A1h	26"W x 20"D	A1s	30"W x 24"D	A1dd
22"W x 19"D	A1i	26"W x 22"D	A1t	30"W x 28"D	A1ee
22"W x 20"D	A1j	26"W x 24"D	A1u	32"W x 26"D	A1ff
22"W x 22"D	A1k	28"W x 18"D	A1v	32"W x 28"D	A1gg

*\*Please note: The polyurethane inner liner is an optional component and not included in the top level cushion assemblies.*




Item	Part No.	Name	Qty
A1a	PUL1709030-2016	LINER,POLYURETHANE,SOL,SOL1,SPCF,SPCG,20"W,16"D	1
A1b	PUL1709030-2018	LINER,POLYURETHANE,SOL,SOL1,SPCF,SPCG,20"W,18"D	1
A1c	PUL1709030-2019	LINER,POLYURETHANE,SOL,SOL1,SPCF,SPCG,20"W,19"D	1
A1d	PUL1709030-2020	LINER,POLYURETHANE,SOL,SOL1,SPCF,SPCG,20"W,20"D	1
A1e	PUL1709030-2022	LINER,POLYURETHANE,SOL,SOL1,SPCF,SPCG,20"W,22"D	1
A1f	PUL1709030-2024	LINER,POLYURETHANE,SOL,SOL1,SPCF,SPCG,20"W,24"D	1
A1g	PUL1709030-2216	LINER,POLYURETHANE,SOL,SOL1,SPCF,SPCG,22"W,16"D	1
A1h	PUL1709030-2218	LINER,POLYURETHANE,SOL,SOL1,SPCF,SPCG,22"W,18"D	1
A1i	PUL1709030-2219	LINER,POLYURETHANE,SOL,SOL1,SPCF,SPCG,22"W,19"D	1
A1j	PUL1709030-2220	LINER,POLYURETHANE,SOL,SOL1,SPCF,SPCG,22"W,20"D	1
A1k	PUL1709030-2222	LINER,POLYURETHANE,SOL,SOL1,SPCF,SPCG,22"W,22"D	1
A1l	PUL1709030-2224	LINER,POLYURETHANE,SOL,SOL1,SPCF,SPCG,22"W,24"D	1
A1m	PUL1709030-2416	LINER,POLYURETHANE,SOL,SOL1,SPCF,SPCG,24"W,16"D	1
A1n	PUL1709030-2418	LINER,POLYURETHANE,SOL,SOL1,SPCF,SPCG,24"W,18"D	1
A1o	PUL1709030-2420	LINER,POLYURETHANE,SOL,SOL1,SPCF,SPCG,24"W,20"D	1
A1p	PUL1709030-2422	LINER,POLYURETHANE,SOL,SOL1,SPCF,SPCG,24"W,22"D	1
A1q	PUL1709030-2424	LINER,POLYURETHANE,SOL,SOL1,SPCF,SPCG,24"W,24"D	1
A1r	PUL1709030-2618	LINER,POLYURETHANE,SOL,SOL1,SPCF,SPCG,26"W,18"D	1
A1s	PUL1709030-2620	LINER,POLYURETHANE,SOL,SOL1,SPCF,SPCG,26"W,20"D	1
A1t	PUL1709030-2622	LINER,POLYURETHANE,SOL,SOL1,SPCF,SPCG,26"W,22"D	1
A1u	PUL1709030-2624	LINER,POLYURETHANE,SOL,SOL1,SPCF,SPCG,26"W,24"D	1
A1v	PUL1709030-2818	LINER,POLYURETHANE,SOL,SOL1,SPCF,SPCG,28"W,18"D	1
A1w	PUL1709030-2820	LINER,POLYURETHANE,SOL,SOL1,SPCF,SPCG,28"W,20"D	1
A1x	PUL1709030-2822	LINER,POLYURETHANE,SOL,SOL1,SPCF,SPCG,28"W,22"D	1
A1y	PUL1709030-2824	LINER,POLYURETHANE,SOL,SOL1,SPCF,SPCG,28"W,24"D	1
A1z	PUL1709030-2828	LINER,POLYURETHANE,SOL,SOL1,SPCF,SPCG,28"W,28"D	1
A1a a	PUL1709030-3018	LINER,POLYURETHANE,SOL,SOL1,SPCF,SPCG,30"W,18"D	1
A1b b	PUL1709030-3020	LINER,POLYURETHANE,SOL,SOL1,SPCF,SPCG,30"W,20"D	1
A1c c	PUL1709030-3022	LINER,POLYURETHANE,SOL,SOL1,SPCF,SPCG,30"W,22"D	1
A1d d	PUL1709030-3024	LINER,POLYURETHANE,SOL,SOL1,SPCF,SPCG,30"W,24"D	1
A1e e	PUL1709030-3028	LINER,POLYURETHANE,SOL,SOL1,SPCF,SPCG,30"W,28"D	1
A1ff	PUL1709030-3226	LINER,POLYURETHANE,SOL,SOL1,SPCF,SPCG,32"W,26"D	1

**Polyurethane Inner Liner - 20"W - 32"W  
Used with Solution, Spectrum Gel,  
and Spectrum Foam Cushions**

Size	Ref#	Size	Ref#	Size	Ref#
20"W x 16"D	A1a	22"W x 24"D	A1l	28"W x 20"D	A1w
20"W x 18"D	A1b	24"W x 16"D	A1m	28"W x 22"D	A1x
20"W x 19"D	A1c	24"W x 18"D	A1n	28"W x 24"D	A1y
20"W x 20"D	A1d	24"W x 20"D	A1o	28"W x 28"D	A1z
20"W x 22"D	A1e	24"W x 22"D	A1p	30"W x 18"D	A1aa
20"W x 24"D	A1f	24"W x 24"D	A1q	30"W x 20"D	A1bb
22"W x 16"D	A1g	26"W x 18"D	A1r	30"W x 22"D	A1cc
22"W x 18"D	A1h	26"W x 20"D	A1s	30"W x 24"D	A1dd
22"W x 19"D	A1i	26"W x 22"D	A1t	30"W x 28"D	A1ee
22"W x 20"D	A1j	26"W x 24"D	A1u	32"W x 26"D	A1ff
22"W x 22"D	A1k	28"W x 18"D	A1v	32"W x 28"D	A1gg

*\*Please note: The polyurethane inner liner is an optional component and not included in the top level cushion assemblies.*






Item	Part No.	Name	Qty
A1gg	PUL1709030-3228	LINER,POLYURETHANE,SOL,SOL1,SPCF,SPCG,32"W,28"D	1

3. Inner Covers

Stealth TRU-Comfort 2 Polyurethane Inner Liner			
Size	Ref#	Size	Ref#
16"W x 14"D	A1a	20"W x 18"D	A1j
16"W x 16"D	A1b	20"W x 20"D	A1k
16"W x 18"D	A1c	20"W x 22"D	A1l
16"W x 20"D	A1d	20"W x 24"D	A1m
18"W x 16"D	A1e	22"W x 18"D	A1n
18"W x 18"D	A1f	22"W x 20"D	A1o
18"W x 20"D	A1g	22"W x 22"D	A1p
18"W x 22"D	A1h	22"W x 24"D	A1q
20"W x 16"D	A1i		



\*Please note: The polyurethane inner liner is an optional component and not included in the top level cushion assemblies.


Item	Part No.	Name	Qty
A1a	PUL1709080-1614	LINER,POLYURETHANE,TC2,ESS,16"W,14"D	1
A1b	PUL1709080-1616	LINER,POLYURETHANE,TC2,ESS,16"W,16"D	1
A1c	PUL1709080-1618	LINER,POLYURETHANE,TC2,ESS,16"W,18"D	1
A1d	PUL1709080-1620	LINER,POLYURETHANE,TC2,ESS,16"W,20"D	1
A1e	PUL1709080-1816	LINER,POLYURETHANE,TC2,ESS,18"W,16"D	1
A1f	PUL1709080-1818	LINER,POLYURETHANE,TC2,ESS,18"W,18"D	1
A1g	PUL1709080-1820	LINER,POLYURETHANE,TC2,ESS,18"W,20"D	1
A1h	PUL1709080-1822	LINER,POLYURETHANE,TC2,ESS,18"W,22"D	1
A1i	PUL1709080-2016	LINER,POLYURETHANE,TC2,ESS,20"W,16"D	1
A1j	PUL1709080-2018	LINER,POLYURETHANE,TC2,ESS,20"W,18"D	1
A1k	PUL1709080-2020	LINER,POLYURETHANE,TC2,ESS,20"W,20"D	1
A1l	PUL1709080-2022	LINER,POLYURETHANE,TC2,ESS,20"W,22"D	1
A1m	PUL1709080-2024	LINER,POLYURETHANE,TC2,ESS,20"W,24"D	1
A1n	PUL1709080-2218	LINER,POLYURETHANE,TC2,ESS,22"W,18"D	1
A1o	PUL1709080-2220	LINER,POLYURETHANE,TC2,ESS,22"W,20"D	1
A1p	PUL1709080-2222	LINER,POLYURETHANE,TC2,ESS,22"W,22"D	1
A1q	PUL1709080-2224	LINER,POLYURETHANE,TC2,ESS,22"W,24"D	1

#### 4. ABS Pans - Solution, Spectrum Gel, Spectrum Foam

**ABS Pan - 10"W - 19"W**  
Used with Stealth Simplicity, Solution, Spectrum Gel, and Spectrum Foam Cushions

Size	Ref#	Size	Ref#	Size	Ref#
10"W x 10"D	A1a	16"W x 18"D	A1j	18"W x 18"D	A1p
10"W x 12"D	A1b	16"W x 19"D	A1k	18"W x 19"D	A1q
12"W x 10"D	A1c	16"W x 20"D	A1k	18"W x 20"D	A1q
12"W x 12"D	A1d	17"W x 16"D	A1l	18"W x 22"D	A1r
14"W x 12"D	A1e	17"W x 17"D	A1l	18"W x 24"D	A1s
14"W x 14"D	A1f	17"W x 18"D	A1m	19"W x 17"D	A1t
14"W x 16"D	A1g	17"W x 19"D	A1m	19"W x 18"D	A1u
14"W x 18"D	A1h	17"W x 20"D	A1n	19"W x 19"D	A1v
16"W x 16"D	A1i	18"W x 16"D	A1o	19"W x 20"D	A1w
16"W x 17"D	A1j	18"W x 17"D	A1o		

\*Sizes shown represent overall cushion size. Actual ABS pan size differs.  
\*\*The ABS pan is an optional component for the cushions listed and is not included with the top level cushion assemblies.




Item	Part No.	Name	Qty
A1a	DWR1010E187A	SOLID SEAT INSERT, 7.5" W X 9.5" D X .125" THK, E-1010-187	1
A1b	DWR1010E187C	SOLID SEAT INSERT, 7.5" W X 11.5" D X .125" THK, E-1010-187	1
A1c	DWR1010E188A	SOLID SEAT INSERT, 9.5" W X 9.5" D X .125" THK, E-1010-188	1
A1d	DWR1010E188C	SOLID SEAT INSERT, 9.5" W X 11.5" D X .125" THK, E-1010-188	1
A1e	DWR1010E189C	SOLID SEAT INSERT, 11.5" W X 11.5" D X .125" THK, E-1010-189	1
A1f	DWR1010E189E	SOLID SEAT INSERT, 11.5" W X 13.5" D X .125" THK, E-1010-189	1
A1g	DWR1010E189G	SOLID SEAT INSERT, 11.5" W X 15.5" D X .125" THK, E-1010-189	1
A1h	DWR1010E189I	SOLID SEAT INSERT, 11.5" W X 17.5" D X .125" THK, E-1010-189	1
A1i	DWR1010E190G	SOLID SEAT INSERT, 13.5" W X 15.5" D X .125" THK, E-1010-190	1
A1j	DWR1010E190I	SOLID SEAT INSERT, 13.5" W X 17.5" D X .125" THK, E-1010-190	1
A1k	DWR1010E190K	SOLID SEAT INSERT, 13.5" W X 19.5" D X .125" THK, E-1010-190	1
A1l	DWR1010E232H	SOLID SEAT INSERT, 14.5" W X 16.5" D X .125" THK, E-1010-232	1
A1m	DWR1010E232J	SOLID SEAT INSERT, 14.5" W X 18.5" D X .125" THK, E-1010-232	1
A1n	DWR1010E190K	SOLID SEAT INSERT, 13.5" W X 19.5" D X .125" THK, E-1010-190	1
A1o	DWR1010E191G	SOLID SEAT INSERT, 15.5" W X 15.5" D X .125" THK, E-1010-191	1
A1p	DWR1010E191I	SOLID SEAT INSERT, 15.5" W X 17.5" D X .125" THK, E-1010-191	1
A1q	DWR1010E191K	SOLID SEAT INSERT, 15.5" W X 19.5" D X .125" THK, E-1010-191	1
A1r	DWR1010E191M	SOLID SEAT INSERT, 15.5" W X 21.5" D X .125" THK, E-1010-191	1
A1s	DWR1010E191O	SOLID SEAT INSERT, 15.5" W X 23.5" D X .125" THK, E-1010-191	1
A1t	DWR1010E234H	SOLID SEAT INSERT, 16.5" W X 16.5" D X .125" THK, E-1010-234	1
A1u	DWR1010E191I	SOLID SEAT INSERT, 15.5" W X 17.5" D X .125" THK, E-1010-191	1
A1v	DWR1010E234J	SOLID SEAT INSERT, 16.5" W X 18.5" D X .125" THK, E-1010-234	1
A1w	DWR1010E192K	SOLID SEAT INSERT, 17.5" W X 19.5" D X .125" THK, E-1010-192	1

#### 4. ABS Pans - Solution, Spectrum Gel, Spectrum Foam

**ABS Pan - 20"W - 32"W**  
Used with Stealth Simplicity, Solution, Spectrum Gel, and Spectrum Foam Cushions

Size	Ref#	Size	Ref#	Size	Ref#
20"W x 16"D	A1a	24"W x 16"D	A1k	28"W x 24"D	A1w
20"W x 18"D	A1b	24"W x 18"D	A1l	28"W x 28"D	A1x
20"W x 19"D	A1c	24"W x 20"D	A1m	30"W x 18"D	A1y
20"W x 20"D	A1c	24"W x 22"D	A1n	30"W x 20"D	A1z
20"W x 22"D	A1d	24"W x 24"D	A1o	30"W x 22"D	A1aa
20"W x 24"D	A1e	26"W x 18"D	A1p	30"W x 24"D	A1bb
22"W x 16"D	A1f	26"W x 20"D	A1q	30"W x 28"D	A1cc
22"W x 18"D	A1g	26"W x 22"D	A1r	32"W x 26"D	A1dd
22"W x 19"D	A1h	26"W x 24"D	A1s	32"W x 28"D	A1ee
22"W x 20"D	A1h	28"W x 18"D	A1t		
22"W x 22"D	A1i	28"W x 20"D	A1u		
22"W x 24"D	A1j	28"W x 22"D	A1v		



\*Sizes shown represent overall cushion size. Actual ABS pan size differs.  
\*\*The ABS pan is an optional component for the cushions listed and is not included with the top level cushion assemblies.



Item	Part No.	Name	Qty
A1a	DWR1010E192G	SOLID SEAT INSERT, 17.5" W X 15.5" D X .125" THK, E-1010-192	1
A1b	DWR1010E192I	SOLID SEAT INSERT, 17.5" W X 17.5" D X .125" THK, E-1010-192	1
A1c	DWR1010E192K	SOLID SEAT INSERT, 17.5" W X 19.5" D X .125" THK, E-1010-192	1
A1d	DWR1010E192M	SOLID SEAT INSERT, 17.5" W X 21.5" D X .125" THK, E-1010-192	1
A1e	DWR1010E192O	SOLID SEAT INSERT, 17.5" W X 23.5" D X .125" THK, E-1010-192	1
A1f	DWR1010E193G	SOLID SEAT INSERT, 19.5" W X 15.5" D X .125" THK, E-1010-193	1
A1g	DWR1010E193I	SOLID SEAT INSERT, 19.5" W X 17.5" D X .125" THK, E-1010-193	1
A1h	DWR1010E193K	SOLID SEAT INSERT, 19.5" W X 19.5" D X .125" THK, E-1010-193	1
A1i	DWR1010E193M	SOLID SEAT INSERT, 19.5" W X 21.5" D X .125" THK, E-1010-193	1
A1j	DWR1010E193O	SOLID SEAT INSERT, 19.5" W X 23.5" D X .125" THK, E-1010-193	1
A1k	DWR1010E194G	SOLID SEAT INSERT, 21.5" W X 15.5" D X .125" THK, E-1010-194	1
A1l	DWR1010E194I	SOLID SEAT INSERT, 21.5" W X 17.5" D X .125" THK, E-1010-194	1
A1m	DWR1010E194K	SOLID SEAT INSERT, 21.5" W X 19.5" D X .125" THK, E-1010-194	1
A1n	DWR1010E194M	SOLID SEAT INSERT, 21.5" W X 21.5" D X .125" THK, E-1010-194	1
A1o	DWR1010E194O	SOLID SEAT INSERT, 21.5" W X 23.5" D X .125" THK, E-1010-194	1
A1p	DWR1010E236I	SOLID SEAT INSERT, 23.5" W X 17.5" D X .125" THK, E-1010-236	1
A1q	DWR1010E236K	SOLID SEAT INSERT, 23.5" W X 19.5" D X .125" THK, E-1010-236	1
A1r	DWR1010E236M	SOLID SEAT INSERT, 23.5" W X 21.5" D X .125" THK, E-1010-236	1
A1s	DWR1010E236O	SOLID SEAT INSERT, 23.5" W X 23.5" D X .125" THK, E-1010-236	1
A1t	DWR1010E240I	SOLID SEAT INSERT, 25.5" W X 17.5" D X .125" THK, E-1010-240	1
A1u	DWR1010E240K	SOLID SEAT INSERT, 25.5" W X 19.5" D X .125" THK, E-1010-240	1
A1v	DWR1010E240M	SOLID SEAT INSERT, 25.5" W X 21.5" D X .125" THK, E-1010-240	1
A1w	DWR1010E240O	SOLID SEAT INSERT, 25.5" W X 23.5" D X .125" THK, E-1010-240	1
A1x	DWR1010E240S	SOLID SEAT INSERT, 25.5" W X 27.5" D X .125" THK, E-1010-240	1
A1y	DWR1010E244I	SOLID SEAT INSERT, 27.5" W X 17.5" D X .125" THK, E-1010-244	1
A1z	DWR1010E244K	SOLID SEAT INSERT, 27.5" W X 19.5" D X .125" THK, E-1010-244	1
A1a	DWR1010E244M	SOLID SEAT INSERT, 27.5" W X 21.5" D X .125" THK, E-1010-244	1
A1b	DWR1010E244O	SOLID SEAT INSERT, 27.5" W X 23.5" D X .125" THK, E-1010-244	1
A1c	DWR1010E244S	SOLID SEAT INSERT, 27.5" W X 27.5" D X .125" THK, E-1010-244	1
A1d	DWR1010E341Q	SOLID SEAT INSERT, 29.5" W X 25.5" D X .125" THK, E-1010-341	1
A1e	DWR1010E341S	SOLID SEAT INSERT, 29.5" W X 27.5" D X .125" THK, E-1010-341	1

5. ACM Board

Complete ACM Board Assy's	
Size	Ref#
16"W x 16"D, 16"W x 18"D	A1
18"W x 16"D, 18"W x 18"D 18"W x 20"D	B1
20"W x 18"D 20"W x 20"D	C1



\*18"W ACM Board Shown

Please note: The ACM boards must be cut to achieve smaller depths shown.






# ACM Boards




Item	Part No.	Name	Qty
A1	CUS142467	ACM BOARD,16"W X 18"D	1
B1	CUS142466	ACM BOARD,18"W X 20"D	1
C1	CUS142465	ACM BOARD, 20"W X 20"D	1

6. Stealth Essence Wedge Assy's



\*Please note: The front half of the ABS pan is glued to the Essence cushion. Wedges fit into the rear and are secured with a hook and loop fastener to the ABS pan. ABS pan non-serviceable.



Item	Part No.	Name	Qty
A1	POP1712340-SM	SEAT WEDGE,POSITIONING KIT,SMALL	1
B1	POP1712340-MD	SEAT WEDGE,POSITIONING KIT,MEDIUM	1
C1	POP1712340-LG	SEAT WEDGE,POSITIONING KIT,LARGE	1